

# Faster. Easier. Reliable. SLE Simplified.

Novum SLE (patent pending) eliminates the questions associated with diatomaceous earth. As a lab manufactured sorbent, you can expect reliable, more consistent results.

## Novum SLE


Reliable synthetic sorbent that excels at batch-to-batch reproducibility due to manufacturing and quality control checks.

## NEW Novum PRO SLE

Same reliable synthetic sorbent with additional clean manufacturing steps to reach low levels of detection for sensitive MS detection, with the same quality reproducibility as Novum SLE for high throughput samples.



# Novum Simplified Liquid Extraction (SLE) Application Guide



- Estrogens from Plasma using Novum PRO
- Testosterone from Plasma using Novum PRO
- Hormone Panel from Plasma using Novum PRO
- Comparison of Synthetic Novum PRO vs Traditional Diatomaceous Earth SLE
- Drugs of Abuse from Urine using Novum
- NSAIDs from Plasma using Novum
- Parabens from Baby Wash using Novum
- COOH-THC from Urine using Novum
- Total Cortisol from Plasma using Novum
- Acrylamide from Coffee
- Automated Clean-up of Cortisone and Prednisolone from Plasma using Novum
- Liquid-Liquid Extraction (LLE) vs. Novum Case Study: Analysis of Monomeric Phenolic Fermentation Inhibitors in Dilute-Acid Plant Hydrolysis

Have Questions?



**chat now**

Quotes, Methods, Tips... We're here to help

[www.phenomenex.com/ChatNow](http://www.phenomenex.com/ChatNow)

Terms and Conditions  
Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Novum is a trademark of Phenomenex.

Novum is patent pending.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.