


# SecurityLINK<sup>SM</sup>

UHPLC Connections in a Click

## SecurityLINK System Connection Points

The SecurityLINK UHPLC fingertight fitting system simplifies your system and column connections while providing consistent performance through torque limiting technology that prevents damage caused by overtightening.



Learn More: [Phenomenex.com/SecurityLINK](http://Phenomenex.com/SecurityLINK)


**phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>

## Installation Instructions

1. Insert SecurityLINK UHPLC fitting into column port



2. Fingertighten until first "CLICK" feedback is received. Connection Complete!



## Ordering Information

### PEEKsil™



PEEKsil Double-Sided 10-32 Fittings with 1/16 in. OD tubing

Part No.	ID (µm)	Length (mm)
AJ1-2111	25	100
AJ1-2121	25	150
AJ1-2141	25	250
AJ1-2151	25	300
AJ1-2171	25	500
AJ1-2191	25	750
AJ1-21A1	25	1000
AJ1-2211	50	100
AJ1-2221	50	150
AJ1-2231	50	200
AJ1-2241	50	250
AJ1-2251	50	300
AJ1-2271	50	500
AJ1-2291	50	750
AJ1-22A1	50	1000
AJ1-2321	75	150
AJ1-2341	75	250
AJ1-2371	75	500
AJ1-23A1	75	1000
AJ1-2411	100	100
AJ1-2421	100	150
AJ1-2441	100	250
AJ1-2471	100	500
AJ1-24A1	100	1000

### Stainless Steel



Stainless Steel Double-Sided 10-32 Fittings with 1/16 in. OD tubing

Part No.	ID (µm)	Length (mm)
AJ1-1421	100	150
AJ1-1441	100	250
AJ1-1461	100	350
AJ1-1471	100	500
AJ1-1481	100	600
AJ1-1521	125	150
AJ1-1541	125	250
AJ1-1561	125	350
AJ1-1571	125	500
AJ1-1581	125	600
AJ1-1621	254	150
AJ1-1641	254	250
AJ1-1661	254	350
AJ1-1671	254	500
AJ1-1681	254	600

#### Terms and Conditions

Subject to Phenomenex standard terms and conditions which can be viewed at [www.phenomenex.com/TermsandConditions](http://www.phenomenex.com/TermsandConditions)

#### Trademarks

SecurityLINK, SecurityCAP, and SecurityGuard are trademarks of Phenomenex. PEEKsil is a trademark of SGE International Pty. Ltd.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2018 Phenomenex, Inc. All rights reserved.

### PEEKsil



PEEKsil Single-Sided 10-32 Fitting with 1/16 in. OD tubing

Part No.	ID (µm)	Length (mm)
AJ1-2224	50	150
AJ1-2274	50	500
AJ1-2294	50	750
AJ1-22A4	50	1000

SecurityLINK tubing material includes a sleeve that provides ID, length and part number information.

**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)