

Zebron Gas Filter Selection Tool

Purity of gas is important for generating high-quality analytical data.

Gas impurities can kill GC columns and detectors and can lead to time intensive troubleshooting.

Zebron Gas Management Filters help protect the GC instrument and extend the column lifetime, minimize stationary phase bleed, ensure low noise, and provide steady, stable baseline.

Select your Gas Management System for GC or GC-MS analysis

What is your
Detector?

FID/FPD/NPD

Fitting


1/4 inch

1/8 inch

Connecting Unit

AG5-2103
4 Position

AG5-2203
4 Position




Your Choice of Zebron Gas Filter




AG5-1030
Hydrocarbon



AG5-1030
Hydrocarbon



AG5-1010
Oxygen



AG5-1020
Moisture



ECD/TCD

1/4 inch


1/8 inch

AG5-2102
2 Position

AG5-2202
2 Position



AG5-1010
Oxygen



AG5-1020
Moisture

MSD

1/4 inch

1/8 inch

AG5-2101
1 Position

AG5-2201
1 Position



AG5-1040
Universal

