



SecurityGuard™ Ultra

Guard Cartridge System for Your Ultra-High Performance Columns

Protect Your Investment

When contaminants and particulates build up at the head of your column, system backpressures can increase dramatically. By simply replacing the SecurityGuard Ultra cartridge, instead of your column, you are able to regain normal operating conditions and reclaim original column performance.

- Extends column lifetimes
- Virtually no change in chromatography
- Easy to use on virtually all manufacturers' UHPLC columns



Accelerated lifetime test using endogenous biological matrix on Kinetex 2.6 μm C18 50 x 4.6 mm ID

* See p. 47-48 for SecurityGuard Ultra ordering information.