

Safety Data Sheet

Date of issue: 19/01/2017 Revision date: : Version: 1.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : 20 ppm SV Test Mix
Product code : AL0-130063
Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Laboratory Use Industrial/Professional use spec : Industrial

For professional use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Phenova

6390 Joyce Dr. Suite 100

80403 Golden, CO - United States T 1-866-942-2978 - F 1-866-283-0269 info@phenova.com - www.phenova.com

1.4. Emergency telephone number

Emergency number : ChemTel Assistance (US/Canada) 1-800-255-3924

ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Flam. Liq. 2 H225 Carc. 2 H351 Aquatic Chronic 3 H412

Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]

Carc.Cat.3; R40 F; R11 R52/53

Full text of R-phrases: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS02



GHS08

Signal word (CLP) : Danger

Hazardous ingredients : Methylene Chloride

Hazard statements (CLP) : H225 - Highly flammable liquid and vapor

H351 - Suspected of causing cancer

H412 - Harmful to aquatic life with long lasting effects

Precautionary statements (CLP) : P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking

19/01/2017 EN (English US) 1/8

Safety Data Sheet

P233 - Keep container tightly closed P273 - Avoid release to the environment

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P303+P361+P353 - IF ON SKIN (or hair): Take off immediately all contaminated clothing.

Rinse skin with water/shower

P308+P313 - IF exposed or concerned: Get medical advice/attention

P403+P235 - Store in a well-ventilated place. Keep cool

No labeling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Methylene Chloride (Component)	(CAS No) 75-09-2 (EC no) 200-838-9 (EC index no) 602-004-00-3	99.98675	Carc. 2, H351
4,4'-DDT (Component)	(CAS No) 50-29-3 (EC no) 200-024-3 (EC index no) 602-045-00-7	0.00265	Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Carc. 2, H351 STOT RE 1, H372 Aquatic Acute 1, H400 (M=100) Aquatic Chronic 1, H410 (M=100)
pyridine (Component) substance with a Community workplace exposure limit	(CAS No) 110-86-1 (EC no) 203-809-9 (EC index no) 613-002-00-7	0.00265	Flam. Liq. 2, H225 Acute Tox. 4 (Oral), H302 Acute Tox. 4 (Dermal), H312 Acute Tox. 4 (Inhalation), H332

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. IF exposed or concerned: Get

medical advice/attention.

First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.

First-aid measures after skin contact : Rinse skin with water/shower. Remove/Take off immediately all contaminated clothing.

First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness

persist.

First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Highly flammable liquid and vapor.

Explosion hazard : May form flammable/explosive vapor-air mixture.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

19/01/2017 EN (English US) 2/8

Safety Data Sheet

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.

Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up in absorbent material. Collect spillage.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Additional hazards when processed : Handle empty containers with care because residual vapors are flammable.

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or

smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. No open flames. No smoking. Use only non-sparking tools. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood.

: Gently wash with plenty of soap and water. Remove/Take off immediately all contaminated

clothing. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

Technical measures : Proper grounding procedures to avoid static electricity should be followed. Ground/bond

container and receiving equipment.

Storage conditions : Keep in fireproof place. Keep container tightly closed. Keep container tightly closed and in a

well-ventilated place. Keep away from any flames or sparking source.

Incompatible materials : Direct sunlight. Heat sources.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

4,4'-DDT (50-29-3)
------------	----------

Hygiene measures

USA OSHA OSHA PEL (TWA) (mg/m³) 1 mg/m³

8.2. Exposure controls

Appropriate engineering controls : Either local exhaust or general room ventilation is usually required.

Personal protective equipment : Avoid all unnecessary exposure. Gloves. Protective clothing. Protective goggles. Safety

glasses.







Hand protection : Wear chemically resistant protective gloves. Wear suitable gloves resistant to chemical

penetration.

Eye protection : Chemical goggles or safety glasses. Safety glasses.

Skin and body protection : Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin

contact.

Respiratory protection : Wear appropriate mask.

Other information : Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid Color Colorless. Odor : characteristic рΗ : No data available Melting point No data available Freezing point No data available Boiling point No data available Flash point : No data available Auto-ignition temperature : No data available

19/01/2017 EN (English US) 3/8

Safety Data Sheet

Decomposition temperature : No data available

Flammability (solid, gas) : Highly flammable liquid and vapor

Relative density : No data available
Solubility : No data available
Explosive properties : No data available
Oxidizing properties : No data available
Explosion limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Highly flammable liquid and vapor. May form flammable/explosive vapor-air mixture.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures. Open flame.

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

May release flammable gases.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

4,4'-DDT (50-29-3)	
LD50 oral rat	87 mg/kg
LD50 dermal rabbit	300 mg/kg
ATE CLP (oral)	87.000 mg/kg body weight
ATE CLP (dermal)	300.000 mg/kg body weight

pyridine (110-86-1)	
LD50 oral rat	> 891 mg/kg (Rat)
LD50 dermal rabbit	1120 mg/kg (Rabbit)
ATE CLP (oral)	500.000 mg/kg body weight
ATE CLP (dermal)	1120.000 mg/kg body weight
ATE CLP (gases)	4500.000 ppmV/4h
ATE CLP (vapors)	11.000 mg/l/4h
ATE CLP (dust, mist)	1.500 mg/l/4h

Methylene Chloride (75-09-2)	
LD50 oral rat	> 2000 mg/kg (Rat; Literature study)
LD50 dermal rabbit	> 2000 mg/kg (Rabbit; Literature study)

Skin corrosion/irritation : Not classified

Based on available data, the classification criteria are not met

Serious eye damage/irritation : Not classified

Based on available data, the classification criteria are not met

Respiratory or skin sensitization : Not classified

Based on available data, the classification criteria are not met

Germ cell mutagenicity : Not classified

Based on available data, the classification criteria are not met

Carcinogenicity : Suspected of causing cancer.

May cause cancer

Reproductive toxicity : Not classified

Based on available data, the classification criteria are not met

Specific target organ toxicity (single exposure) : Not classified

Based on available data, the classification criteria are not met

19/01/2017 EN (English US) 4/8

Safety Data Sheet

Specific target organ toxicity (repeated : Not classified

exposure) Based on available data, the classification criteria are not met

Aspiration hazard : Not classified

Based on available data, the classification criteria are not met

Potential Adverse human health effects and

symptoms

: Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1. Toxicity

Ecology - water : Harmful to aquatic life with long lasting effects.

4,4'-DDT (50-29-3)		
LC50 fish 1	0.01 mg/l Pimephales promelas (fathead minnow) 96 h	
LC50 other aquatic organisms 1	0.0034 mg/l Oncorhynchus mykiss (rainbow trouit) 96 h	
EC50 Daphnia 1	0.00108 mg/l Immolbilization - Daphnia magna (Water flea) 48 h	
LC50 fish 2	0.01 mg/l Lepomis macrochirus (Bluegill) 96 h	
LOEC (acute)	150 mg/l Oncorhynchus mykiss (rainbow trout) 3 d	
NOEC (acute)	113 mg/l Oncorhynchus mykiss (rainbow trout) 3 d	
pyridine (110-86-1)		
LC50 fish 1	4.6 mg/l (LC50; 96 h)	
EC50 Daphnia 2	495 mg/l (EC50; 48 h)	
Methylene Chloride (75-09-2)		
LC50 fish 1	193 mg/l (LC50; 96 h; Pimephales promelas)	
EC50 Daphnia 1	168.2 mg/l (EC50; 48 h)	

12.2. Persistence and degradability

20 ppm SV Test Mix		
Persistence and degradability	May cause long-term adverse effects in the environment.	
pyridine (110-86-1)		
Persistence and degradability	Readily biodegradable in water. Non degradable in the soil. Biodegradable in the soil under anaerobic conditions.	
Biochemical oxygen demand (BOD)	1.15 g O /g substance	
Chemical oxygen demand (COD)	0.05 g O /g substance	
ThOD	2.23 g O /g substance	
BOD (% of ThOD)	0.52	
Methylene Chloride (75-09-2)	·	

Wethylene Chloride (75-05-2)

Persistence and degradability Not readily biodegradable in water. Biodegradable in the soil.

12.3. Bioaccumulative potential

20 ppm SV Test Mix		
Not established.		
4,4'-DDT (50-29-3)		
46670 Oncorhynchus mykiss (rainbow trout) 20 d		
6.91		
pyridine (110-86-1)		
0.65 - 1.04 (Experimental value)		
Low potential for bioaccumulation (Log Kow < 4).		
Methylene Chloride (75-09-2)		
2 - 40 (BCF)		
1.25 (Experimental value)		

Bioaccumulative potential 12.4. Mobility in soil

pyridine (110-86-1)		
Surface tension	0.038 N/m (20 °C)	
Methylene Chloride (75-09-2)		
Surface tension	0.028 N/m (20 °C)	
Ecology - soil	May be harmful to plant growth, blooming and fruit formation.	

Low potential for bioaccumulation (BCF < 500).

12.5. Results of PBT and vPvB assessment

No additional information available

19/01/2017 EN (English US) 5/8

Safety Data Sheet

12.6. Other adverse effects

Additional information : Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

Additional information : Handle empty containers with care because residual vapors are flammable.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

 UN-No. (ADR)
 : 1993

 UN-No. (IATA)
 : 1993

 UN-No. (IMDG)
 : 1993

 UN-No.(ADN)
 : 1993

14.2. UN proper shipping name

Proper Shipping Name (ADR) : FLAMMABLE LIQUID, N.O.S.
Proper Shipping Name (IATA) : Flammable liquid, n.o.s.
Proper Shipping Name (IMDG) : FLAMMABLE LIQUID, N.O.S.
Proper Shipping Name (ADN) : FLAMMABLE LIQUID, N.O.S.

Transport document description (ADR) : UN 1993 FLAMMABLE LIQUID, N.O.S., 3, II, (D/E)

14.3. Packing group

 Class (ADR)
 : 3

 Classification code (ADR)
 : F1

 Class (IATA)
 : 3

 Class (IMDG)
 : 3

 Class (ADN)
 : 3

 Classification code (ADN)
 : F1

 Hazard labels (ADR)
 : 3



Hazard labels (IATA) : 3



Hazard labels (IMDG) : 3



Hazard labels (ADN) : 3

14.4. Packing group	o
---------------------	---

 Packing group (ADR)
 : II

 Packing group (IATA)
 : II

 Packing group (IMDG)
 : II

 Packing group (ADN)
 : II

14.5. Environmental hazards

Other information : No supplementary information available.

14.6. Special precautions for user

14.6.1. Overland transport

Hazard identification number (Kemler No.) : 33

19/01/2017 EN (English US) 6/8

Safety Data Sheet

Classification code (ADR) : F1

Orange plates :

33 1993

Special provision (ADR) : 274, 601, 640D

Transport category (ADR) : 2
Tunnel restriction code (ADR) : D/E
Limited quantities (ADR) : 11
Excepted quantities (ADR) : E2

14.6.2. Transport by sea

Special provision (IMDG) : 274
Limited quantities (IMDG) : 1 L
Excepted quantities (IMDG) : E2
Packing instructions (IMDG) : P001
IBC packing instructions (IMDG) : IBC02
Tank instructions (IMDG) : T7

Tank special provisions (IMDG) : TP1, TP8, TP28

 EmS-No. (Fire)
 : F-E

 EmS-No. (Spillage)
 : S-E

 Stowage category (IMDG)
 : B

14.6.3. Air transport

CAO packing instructions (IATA) : 364 CAO max net quantity (IATA) : 60L PCA packing instructions (IATA) : 353 PCA Limited quantities (IATA) : Y341 PCA limited quantity max net quantity (IATA) : 1L PCA max net quantity (IATA) : 5L PCA Excepted quantities (IATA) : E2 Special provision (IATA) : A3 ERG code (IATA) : 3H

14.6.4. Inland waterway transport

Special provision (ADN) : 274, 601, 640D

Limited quantities (ADN) : 1 L

Excepted quantities (ADN) : E2

Carriage permitted (ADN) : T

Equipment required (ADN) : PP, EX, A

Ventilation (ADN) : VE01

Number of blue cones/lights (ADN) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

19/01/2017 EN (English US) 7/8

Safety Data Sheet

REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and Data sources

mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Other information : None.

PHV SDS EU

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.

19/01/2017 EN (English US) 8/8