


96-Well Collection Plates from Phenomenex

- Available in a variety of options including round well, square well, and multiple volumes
- Perfect for general lab use and SPE procedures
- Footprints conform to ANSI/SBS standards for easy robotic programming



ROUND WELL

Part Number: **AH0-9341**
Well Volume: 0.5 mL
Well Shape: Round
Well Bottom: Conical
Material: Polypropylene


Part Number: **AH0-9333**
Well Volume: 0.5 mL
Well Shape: Round
Well Bottom: Conical
Material: Polypropylene
Special Notes: **Part No. AH0-9333 is Sterile**

Recommended Sealing Mats

AH0-8631 Pierceable, 96-Round Well 7 mm, Silicone

AH0-8632 Pre-Slit, 96-Round Well 7 mm, Silicone

AH0-7362 Sealing Tape Pad



All dimensions are measured in mm.

ROUND WELL


Part Number: **AH0-7279**
Well Volume: 1 mL
Well Shape: Round
Well Bottom: Round
Material: Polypropylene

Recommended Sealing Mats

AH0-8631 Pierceable, 96-Round Well 7 mm, Silicone

AH0-8632 Pre-Slit, 96-Round Well 7 mm, Silicone

AH0-7362 Sealing Tape Pad



All dimensions are measured in mm.

ROUND WELL

Part Number: **AH0-9332**
Well Volume: 1.2 mL
Well Shape: Round
Well Bottom: Conical
Material: Polypropylene


Part Number: **AH0-9333**
Well Volume: 0.5 mL
Well Shape: Round
Well Bottom: Conical
Material: Polypropylene
Special Notes: **Part No. AH0-9333 is Sterile**

Recommended Sealing Mats

AH0-8633 Pierceable, 96-Round Well 8 mm, Silicone

AH0-8634 Pre-Slit 96-Round Well 8 mm, Silicone

AH0-7362 Sealing Tape Pad



All dimensions are measured in mm.

ROUND WELL


Part Number: **AH0-8636**
Well Volume: 2 mL
Well Shape: Round
Well Bottom: Round
Material: Polypropylene

Recommended Sealing Mats

AH0-8631 Pierceable, 96-Round Well 8 mm, Silicone

AH0-8632 Pre-Slit 96-Round Well 8 mm, Silicone

AH0-7362 Sealing Tape Pad




All dimensions are measured in mm.

SQUARE WELL

Part Number: **AH0-7192**
Well Volume: 350 μ L
Well Shape: Square
Well Bottom: Conical
Material: Polypropylene

Recommended Sealing Mats

AH0-8597 Pierceable, 96-Square Well, Silicone
AH0-8598 Pre-Slit, 96-Square Well, Silicone
AH0-8199 Pierceable, 96-Square Well, Santoprene™
AH0-7195 Pierceable, 96-Square Well, Ethylene Vinyl Acetate (EVA)
AH0-7362 Sealing Tape Pad




All dimensions are measured in mm.

SQUARE WELL

Part Number: **AH0-7193**
Well Volume: 1 mL
Well Shape: Square
Well Bottom: Conical
Material: Polypropylene

Recommended Sealing Mats

AH0-8597 Pierceable, 96-Square Well, Silicone
AH0-8598 Pre-Slit, 96-Square Well, Silicone
AH0-8199 Pierceable, 96-Square Well, Santoprene™
AH0-7195 Pierceable, 96-Square Well, Ethylene Vinyl Acetate (EVA)
AH0-7362 Sealing Tape Pad




All dimensions are measured in mm.

SQUARE WELL

Part Number: **AH0-7194**
Well Volume: 2 mL
Well Shape: Square
Well Bottom: Conical
Material: Polypropylene

Recommended Sealing Mats

AH0-8597 Pierceable, 96-Square Well, Silicone
AH0-8598 Pre-Slit, 96-Square Well, Silicone
AH0-8199 Pierceable, 96-Square Well, Santoprene™
AH0-7195 Pierceable, 96-Square Well, Ethylene Vinyl Acetate (EVA)
AH0-7362 Sealing Tape Pad



Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thainfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☉ **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
www.phenomenex.com/chat



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Trademarks

Santoprene is a trademark of Exxon Mobil Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.