

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Date of issue: 08/07/2019

Revision date: 08/07/2019

Version: 1.0

SECTION 1: Identification

1.1. Identification

Product form : Mixture
 Product name : Custom 8270 Short List
 Product code : AL0-130837

1.2. Recommended use and restrictions on use

No additional information available

1.3. Supplier

Phenova
 6390 Joyce Dr. Suite 100
 Golden, CO 80403 - United States
 T 1-866-942-2978 - F 1-866-283-0269
info@phenova.com - www.phenova.com

1.4. Emergency telephone number

Emergency number : ChemTel Assistance (US/Canada) 1-800-255-3924
 ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazard(s) identification

2.1. Classification of the substance or mixture

GHS US classification

Flammable liquids Category 4	H227	Combustible liquid
Skin sensitization, Category 1	H317	May cause an allergic skin reaction
Germ cell mutagenicity Category 1B	H340	May cause genetic defects
Carcinogenicity Category 1B	H350	May cause cancer
Reproductive toxicity Category 1B	H360	May damage fertility or the unborn child

Full text of H statements : see section 16

2.2. GHS Label elements, including precautionary statements

GHS US labeling

Hazard pictograms (GHS-US) :



Signal word (GHS-US) :

Danger

Hazard statements (GHS-US) :

H227 - Combustible liquid
 H317 - May cause an allergic skin reaction
 H340 - May cause genetic defects
 H350 - May cause cancer
 H360 - May damage fertility or the unborn child

Precautionary statements (GHS-US) :

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
 P261 - Avoid breathing dust/fume/gas/mist/vapors/spray.
 P272 - Contaminated work clothing must not be allowed out of the workplace
 P280 - Wear protective gloves/protective clothing/eye protection/face protection.
 P308+P313 - If exposed or concerned: Get medical advice/attention.
 P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.
 P363 - Wash contaminated clothing before reuse.
 P370+P378 - In case of fire: Use media other than water to extinguish.
 P403+P235 - Store in a well-ventilated place. Keep cool.
 P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

2.3. Other hazards which do not result in classification

No additional information available

2.4. Unknown acute toxicity (GHS US)

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	Conc.
Methylene Chloride (Component)	(CAS-No.) 75-09-2	95.2
atrazine (Component)	(CAS-No.) 1912-24-9	0.2
Benzo(b)fluoranthene (Component)	(CAS-No.) 205-99-2	0.2
benzo[a]anthracene (Component)	(CAS-No.) 56-55-3	0.2
Bis(2-ethylhexyl) phthalate (Component)	(CAS-No.) 117-81-7	0.2
chrysene (Component)	(CAS-No.) 218-01-9	0.2
benzo[k]fluoranthene (Component)	(CAS-No.) 207-08-9	0.2
benzo[a]pyrene (Component)	(CAS-No.) 50-32-8	0.2
dibenz(a,h)anthracene (Component)	(CAS-No.) 53-70-3	0.2
hexachlorobuta-1,3-diene (Component)	(CAS-No.) 87-68-3	0.2
naphthalene (Component)	(CAS-No.) 91-20-3	0.2
nitrobenzene (Component)	(CAS-No.) 98-95-3	0.2
hexachloroethane (Component)	(CAS-No.) 67-72-1	0.2
2,3,4,5,6-pentachlorophenol (Component)	(CAS-No.) 87-86-5	0.2
hexachlorobenzene (Component)	(CAS-No.) 118-74-1	0.2
indeno(1,2,3-cd)pyrene (Component)	(CAS-No.) 193-39-5	0.2
N-Nitrosodimethylamine (Component)	(CAS-No.) 62-75-9	0.2

Full text of hazard classes and H-statements : see section 16

SECTION 4: First-aid measures

4.1. Description of first aid measures

First-aid measures general	: Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Allow victim to breathe fresh air. Allow the victim to rest.
First-aid measures after skin contact	: Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse.
First-aid measures after eye contact	: Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persists.
First-aid measures after ingestion	: Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects (acute and delayed)

Potential Adverse human health effects and symptoms	: Based on available data, the classification criteria are not met.
Symptoms/effects	: Not expected to present a significant hazard under anticipated conditions of normal use.

4.3. Immediate medical attention and special treatment, if necessary

No additional information available

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

SECTION 5: Fire-fighting measures

5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.
Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Specific hazards arising from the chemical

No additional information available

5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.
Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect spillage. Store away from other materials.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep only in the original container in a cool, well ventilated place away from : Keep container closed when not in use.
Incompatible products : Strong bases. Strong acids.
Incompatible materials : Sources of ignition. Direct sunlight.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Custom 8270 Short List		
ACGIH	Local name	Dichloromethane
ACGIH	ACGIH TWA (ppm)	50 ppm
ACGIH	Remark (ACGIH)	COHb-emia; CNS impair
ACGIH	Regulatory reference	ACGIH 2018
OSHA	Remark (OSHA)	(2) See Table Z-2.
OSHA	Regulatory reference (US-OSHA)	OSHA
atrazine (1912-24-9)		
ACGIH	Local name	Atrazine
ACGIH	ACGIH TWA (mg/m ³)	2 mg/m ³ (Inhalable fraction)

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

atrazine (1912-24-9)		
ACGIH	Remark (ACGIH)	CNS convul; A3 (Confirmed Animal Carcinogen with Unknown Relevance to Humans: The agent is carcinogenic in experimental animals at a relatively high dose, by route(s) of administration, at site(s), of histologic type(s), or by mechanism(s) that may not be relevant to worker exposure. Available epidemiologic studies do not confirm an increased risk of cancer in exposed humans. Available evidence does not suggest that the agent is likely to cause cancer in humans except under uncommon or unlikely routes or levels of exposure)
ACGIH	Regulatory reference	ACGIH 2018
benzo[a]anthracene (56-55-3)		
Not applicable		
Benzo(b)fluoranthene (205-99-2)		
Not applicable		
benzo[k]fluoranthene (207-08-9)		
Not applicable		
benzo[a]pyrene (50-32-8)		
Not applicable		
chrysene (218-01-9)		
Not applicable		
dibenz(a,h)anthracene (53-70-3)		
Not applicable		
Bis(2-ethylhexyl) phthalate (117-81-7)		
ACGIH	Local name	Di(2-ethylhexyl)phthalate (DEHP)
ACGIH	ACGIH TWA (mg/m ³)	5 mg/m ³
ACGIH	Remark (ACGIH)	LRT irr
ACGIH	Regulatory reference	ACGIH 2018
OSHA	OSHA PEL (TWA) (mg/m ³)	5 mg/m ³
OSHA	Regulatory reference (US-OSHA)	OSHA
hexachlorobenzene (118-74-1)		
ACGIH	Local name	Hexachlorobenzene
ACGIH	ACGIH TWA (mg/m ³)	0.002 mg/m ³
ACGIH	Remark (ACGIH)	Porphyrin eff; Skin dam; CNS impair
ACGIH	Regulatory reference	ACGIH 2018
hexachlorobuta-1,3-diene (87-68-3)		
ACGIH	Local name	Hexachlorobutadiene
ACGIH	ACGIH TWA (ppm)	0.02 ppm
ACGIH	Remark (ACGIH)	Kidney dam
ACGIH	Regulatory reference	ACGIH 2018
hexachloroethane (67-72-1)		
ACGIH	Local name	Hexachloroethane
ACGIH	ACGIH TWA (ppm)	1 ppm
ACGIH	Remark (ACGIH)	Liver & kidney dam
ACGIH	Regulatory reference	ACGIH 2018
OSHA	OSHA PEL (TWA) (mg/m ³)	10 mg/m ³
OSHA	OSHA PEL (TWA) (ppm)	1 ppm
OSHA	Regulatory reference (US-OSHA)	OSHA

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

indeno(1,2,3-cd)pyrene (193-39-5)		
Not applicable		
naphthalene (91-20-3)		
ACGIH	Local name	Naphthalene
ACGIH	ACGIH TWA (ppm)	10 ppm (Naphthalene; USA; Time-weighted average exposure limit 8 h; TLV - Adopted Value)
ACGIH	Remark (ACGIH)	Hematologic eff; URT & eye irr; Skin; A3 (Confirmed Animal Carcinogen with Unknown Relevance to Humans: The agent is carcinogenic in experimental animals at a relatively high dose, by route(s) of administration, at site(s), of histologic type(s), or by mechanism(s) that may not be relevant to worker exposure. Available epidemiologic studies do not confirm an increased risk of cancer in exposed humans. Available evidence does not suggest that the agent is likely to cause cancer in humans except under uncommon or unlikely routes or levels of exposure)
ACGIH	Regulatory reference	ACGIH 2018
OSHA	OSHA PEL (TWA) (mg/m ³)	50 mg/m ³
OSHA	OSHA PEL (TWA) (ppm)	10 ppm
OSHA	Regulatory reference (US-OSHA)	OSHA
nitrobenzene (98-95-3)		
ACGIH	Local name	Nitrobenzene
ACGIH	ACGIH TWA (ppm)	1 ppm (Nitrobenzene; USA; Time-weighted average exposure limit 8 h; TLV - Adopted Value)
ACGIH	Remark (ACGIH)	MeHb-emia
ACGIH	Regulatory reference	ACGIH 2018
OSHA	OSHA PEL (TWA) (mg/m ³)	5 mg/m ³
OSHA	OSHA PEL (TWA) (ppm)	1 ppm
OSHA	Regulatory reference (US-OSHA)	OSHA
N-Nitrosodimethylamine (62-75-9)		
Not applicable		
2,3,4,5,6-pentachlorophenol (87-86-5)		
ACGIH	Local name	Pentachlorophenol
ACGIH	ACGIH TWA (mg/m ³)	0.5 mg/m ³ (Inhalable fraction and vapor)
ACGIH	ACGIH STEL (mg/m ³)	1 mg/m ³ (Inhalable fraction and vapor)
ACGIH	Remark (ACGIH)	URT & eye irr; CNS & card impair; Skin; A3 (Confirmed Animal Carcinogen with Unknown Relevance to Humans: The agent is carcinogenic in experimental animals at a relatively high dose, by route(s) of administration, at site(s), of histologic type(s), or by mechanism(s) that may not be relevant to worker exposure. Available epidemiologic studies do not confirm an increased risk of cancer in exposed humans. Available evidence does not suggest that the agent is likely to cause cancer in humans except under uncommon or unlikely routes or levels of exposure); BEI
ACGIH	Regulatory reference	ACGIH 2018
OSHA	OSHA PEL (TWA) (mg/m ³)	0.5 mg/m ³
OSHA	Regulatory reference (US-OSHA)	OSHA
Methylene Chloride (75-09-2)		
ACGIH	Local name	Dichloromethane
ACGIH	ACGIH TWA (ppm)	50 ppm

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Methylene Chloride (75-09-2)		
ACGIH	Remark (ACGIH)	COHb-emia; CNS impair
ACGIH	Regulatory reference	ACGIH 2018
OSHA	Remark (OSHA)	(2) See Table Z-2.
OSHA	Regulatory reference (US-OSHA)	OSHA

8.2. Appropriate engineering controls

No additional information available

8.3. Individual protection measures/Personal protective equipment

Personal protective equipment:

Avoid all unnecessary exposure.

Hand protection:

Wear protective gloves.

Eye protection:

Chemical goggles or safety glasses

Respiratory protection:

Wear appropriate mask

Other information:

Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
	: Colorless
	: characteristic
Odor threshold	: No data available
pH	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Non flammable.
Vapor pressure	: No data available
Relative vapor density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Log Pow	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: No data available
Oxidizing properties	: No data available

9.2. Other information

No additional information available

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Not established.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

fume. Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Bis(2-ethylhexyl) phthalate (117-81-7)

LD50 dermal rabbit	19800 mg/kg body weight (24 h, Rabbit, Experimental value, Dermal)
ATE US (dermal)	19800 mg/kg body weight

hexachlorobenzene (118-74-1)

LD50 oral rat	10000 mg/kg (Rat, Oral)
ATE US (oral)	10000 mg/kg body weight

hexachlorobuta-1,3-diene (87-68-3)

LD50 oral rat	90 mg/kg (Rat, Oral)
LD50 dermal rabbit	1211 mg/kg (Rabbit, Dermal)
ATE US (oral)	90 mg/kg body weight
ATE US (dermal)	1211 mg/kg body weight

hexachloroethane (67-72-1)

LD50 oral rat	4460 mg/kg (Rat, Oral)
LD50 dermal rabbit	32000 mg/kg (Rabbit, Dermal)
ATE US (oral)	4460 mg/kg body weight
ATE US (dermal)	32000 mg/kg body weight

naphthalene (91-20-3)

LD50 oral rat	> 1100 mg/kg (Rat)
LD50 dermal rat	> 2500 mg/kg (Rat)
LD50 dermal rabbit	> 20000 mg/kg (Rabbit)
ATE US (oral)	500 mg/kg body weight

nitrobenzene (98-95-3)

LD50 oral rat	640 mg/kg (Rat; Experimental value; 588 mg/kg bodyweight; Rat)
LD50 dermal rabbit	760 mg/kg body weight (Rabbit; Experimental value)
ATE US (oral)	100 mg/kg body weight
ATE US (dermal)	760 mg/kg body weight
ATE US (gases)	700 ppmV/4h
ATE US (vapors)	3 mg/l/4h
ATE US (dust, mist)	0.5 mg/l/4h

N-Nitrosodimethylamine (62-75-9)

LD50 oral rat	37 mg/kg (Rat)
LC50 inhalation rat (mg/l)	0.24 mg/l/4h (Rat)
LC50 inhalation rat (ppm)	78 ppm/4h (Rat)

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

N-Nitrosodimethylamine (62-75-9)	
ATE US (oral)	37 mg/kg body weight
ATE US (gases)	78 ppmV/4h
ATE US (vapors)	0.24 mg/l/4h
ATE US (dust, mist)	0.24 mg/l/4h
2,3,4,5,6-pentachlorophenol (87-86-5)	
LC50 inhalation rat (mg/l)	355 mg/m ³ (Rat, Literature, Inhalation)
ATE US (oral)	100 mg/kg body weight
ATE US (dermal)	300 mg/kg body weight
ATE US (gases)	100 ppmV/4h
ATE US (vapors)	0.5 mg/l/4h
ATE US (dust, mist)	0.05 mg/l/4h
Methylene Chloride (75-09-2)	
LD50 oral rat	> 2000 mg/kg body weight (OECD 401: Acute Oral Toxicity, Rat, Male / female, Experimental value, Oral)
LD50 dermal rat	> 2000 mg/kg body weight (OECD 402: Acute Dermal Toxicity, 24 h, Rat, Male / female, Experimental value, Dermal)
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitization	: May cause an allergic skin reaction.
Germ cell mutagenicity	: May cause genetic defects.
Carcinogenicity	: May cause cancer.
benzo[a]anthracene (56-55-3)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
Benzo(b)fluoranthene (205-99-2)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
benzo[k]fluoranthene (207-08-9)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
benzo[a]pyrene (50-32-8)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
dibenz(a,h)anthracene (53-70-3)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
Bis(2-ethylhexyl) phthalate (117-81-7)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
hexachlorobenzene (118-74-1)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
hexachloroethane (67-72-1)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
indeno(1,2,3-cd)pyrene (193-39-5)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
naphthalene (91-20-3)	
IARC group	2B - Possibly carcinogenic to humans
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
nitrobenzene (98-95-3)	
IARC group	2B - Possibly carcinogenic to humans
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
N-Nitrosodimethylamine (62-75-9)	
IARC group	2A - Probably carcinogenic to humans
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen
2,3,4,5,6-pentachlorophenol (87-86-5)	
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Methylene Chloride (75-09-2)	
IARC group	2A - Probably carcinogenic to humans
National Toxicology Program (NTP) Status	Reasonably anticipated to be Human Carcinogen

Reproductive toxicity : May damage fertility or the unborn child.
 Specific target organ toxicity – single exposure : Not classified

Specific target organ toxicity – repeated exposure : Not classified

Aspiration hazard : Not classified

Potential Adverse human health effects and symptoms : Based on available data, the classification criteria are not met.

Symptoms/effects : Not expected to present a significant hazard under anticipated conditions of normal use.

SECTION 12: Ecological information

12.1. Toxicity

benzo[a]anthracene (56-55-3)	
LC50 fish 1	0.0018 mg/l (65 h, Pimephales promelas, Lethal)
EC50 Daphnia 1	0.01 mg/l (96 h, Daphnia pulex, Static system)

benzo[a]pyrene (50-32-8)	
LC50 fish 1	0.0056 mg/l (38 h, Pimephales promelas, Lethal)

hexachlorobenzene (118-74-1)	
LC50 fish 1	2.3 mg/l (96 h, Salmo gairdneri)
EC50 Daphnia 1	> 0.03 mg/l (24 h, Daphnia magna)

hexachlorobuta-1,3-diene (87-68-3)	
LC50 fish 1	0.25 mg/l (96 h, Salmo gairdneri)
EC50 other aquatic organisms 1	0.21 mg/l (96 h, Lymnaea sp.)

hexachloroethane (67-72-1)	
LC50 fish 1	0.84 mg/l (96 h, Salmo gairdneri)
EC50 Daphnia 1	1.4 mg/l (Daphnia magna)

naphthalene (91-20-3)	
EC50 Daphnia 1	2.16 mg/l (EC50; 48 h; Daphnia magna)
LC50 fish 2	0.11 mg/l (LC50; 96 h; Oncorhynchus mykiss)
Threshold limit algae 1	0.4 mg/l (EC50; 72 h; Skeletonema costatum)

nitrobenzene (98-95-3)	
LC50 fish 1	4.3 mg/l (LC50; OECD 203: Fish, Acute Toxicity Test; 48 h; Oryzias latipes)
EC50 Daphnia 1	35 mg/l (Other, 48 h, Daphnia magna, Static system, Fresh water, Experimental value, Locomotor effect)

2,3,4,5,6-pentachlorophenol (87-86-5)	
LC50 fish 1	0.052 mg/l (96 h, Salmo gairdneri)
EC50 Daphnia 1	0.01 - 0.36 mg/l (48 h, Daphnia magna)

Methylene Chloride (75-09-2)	
LC50 fish 1	193 mg/l (96 h, Pimephales promelas, Flow-through system, Fresh water, Experimental value)
EC50 Daphnia 1	168.2 mg/l (48 h, Daphnia magna)

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

12.2. Persistence and degradability

Custom 8270 Short List	
Persistence and degradability	Not established.
atrazine (1912-24-9)	
Persistence and degradability	Biodegradability in soil: no data available. Not readily biodegradable in water.
benzo[a]anthracene (56-55-3)	
Persistence and degradability	Biodegradability in soil: no data available. Inhibits biodegradation processes in the soil. Not readily biodegradable in water.
ThOD	2.95 g O ₂ /g substance
Benzo(b)fluoranthene (205-99-2)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
ThOD	2.92 g O ₂ /g substance
benzo[k]fluoranthene (207-08-9)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
ThOD	2.92 g O ₂ /g substance
benzo[a]pyrene (50-32-8)	
Persistence and degradability	Biodegradable in the soil. Not readily biodegradable in water.
Chemical oxygen demand (COD)	2.92 g O ₂ /g substance
ThOD	2.92 g O ₂ /g substance
chrysene (218-01-9)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
dibenz(a,h)anthracene (53-70-3)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
Bis(2-ethylhexyl) phthalate (117-81-7)	
Persistence and degradability	Biodegradable in the soil. Readily biodegradable in water.
hexachlorobenzene (118-74-1)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water. Not easily biodegradable in water in anaerobic conditions.
hexachlorobuta-1,3-diene (87-68-3)	
Persistence and degradability	Biodegradability in soil: no data available. Readily biodegradable in water.
hexachloroethane (67-72-1)	
Persistence and degradability	Readily biodegradable in water.
indeno(1,2,3-cd)pyrene (193-39-5)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
ThOD	2.9 g O ₂ /g substance
naphthalene (91-20-3)	
Persistence and degradability	Readily biodegradable in water. Forming sediments in water. Biodegradable in the soil. Adsorbs into the soil. Photolysis in the air.
Biochemical oxygen demand (BOD)	0 g O ₂ /g substance
Chemical oxygen demand (COD)	0.22 g O ₂ /g substance
ThOD	2.99 g O ₂ /g substance
nitrobenzene (98-95-3)	
Persistence and degradability	Not readily biodegradable in water. Biodegradable in the soil. Low potential for adsorption in soil.
Biochemical oxygen demand (BOD)	0 g O ₂ /g substance
ThOD	1.95 g O ₂ /g substance
BOD (% of ThOD)	0

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

N-Nitrosodimethylamine (62-75-9)	
Persistence and degradability	Not readily biodegradable in water. Photolysis in water. Photolysis in the air.
2,3,4,5,6-pentachlorophenol (87-86-5)	
Persistence and degradability	Non degradable in the soil. Not readily biodegradable in water.
Methylene Chloride (75-09-2)	
Persistence and degradability	Biodegradable in the soil. Not readily biodegradable in water.
12.3. Bioaccumulative potential	
Custom 8270 Short List	
Bioaccumulative potential	Not established.
atrazine (1912-24-9)	
BCF fish 1	3 - 4 (Cyprinus carpio)
BCF fish 2	3 - 10 (Pisces)
BCF other aquatic organisms 1	52 (24 h, Chlorella sp.)
BCF other aquatic organisms 2	10 - 83 (Algae)
Log Pow	2.64
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).
benzo[a]anthracene (56-55-3)	
BCF fish 1	350 (72 h, Leuciscus idus)
BCF other aquatic organisms 1	1106 (24 h, Daphnia pulex)
BCF other aquatic organisms 2	18000 (192 h, Crassostrea sp.)
Log Pow	5.61 - 5.79
Bioaccumulative potential	High potential for bioaccumulation (BCF > 5000).
Benzo(b)fluoranthene (205-99-2)	
BCF other aquatic organisms 1	2800 (168 h, Lamellibranchiata)
Log Pow	6.57
Bioaccumulative potential	High potential for bioaccumulation (Log Kow > 5).
benzo[k]fluoranthene (207-08-9)	
BCF fish 1	8750 (Pisces, QSAR)
BCF other aquatic organisms 1	0.0013 mg/kg (Algae, Dry weight)
BCF other aquatic organisms 2	37000 (Mytilus edulis)
Log Pow	6.84
Bioaccumulative potential	High potential for bioaccumulation (BCF > 5000).
benzo[a]pyrene (50-32-8)	
BCF fish 1	480 (72 h, Leuciscus idus)
BCF fish 2	70.7 (168 h, Salmo salar, Eggs)
BCF other aquatic organisms 1	3000 (192 h, Crassostrea sp.)
BCF other aquatic organisms 2	1.5 (24 h, Daphnia magna)
Log Pow	5.97 - 6.06
Bioaccumulative potential	High potential for bioaccumulation (Log Kow > 5).
chrysene (218-01-9)	
BCF other aquatic organisms 1	4440 (180 day(s), Lamellibranchiata, Literature study, Chronic)
Log Pow	5.81 - 5.86 (Experimental value)
Bioaccumulative potential	High potential for bioaccumulation (Log Kow > 5).
dibenz(a,h)anthracene (53-70-3)	
Log Pow	5.97 - 6.84
Bis(2-ethylhexyl) phthalate (117-81-7)	
BCF fish 1	155 - 886 (56 day(s), Pimephales promelas, Literature study)
Log Pow	7.68 (Experimental value, Other)
Bioaccumulative potential	High potential for bioaccumulation (Log Kow > 5).
hexachlorobenzene (118-74-1)	
BCF fish 1	20000 (Salmo gairdneri, Test duration: 8 weeks)
BCF fish 2	30000 (Cyprinus carpio, Test duration: 8 weeks)

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

hexachlorobenzene (118-74-1)	
BCF other aquatic organisms 1	25000 (Algae)
BCF other aquatic organisms 2	1130 (720 h, Daphnia magna)
Log Pow	5.73 - 6.39 (Experimental value)
Bioaccumulative potential	High potential for bioaccumulation (BCF > 5000).
hexachlorobuta-1,3-diene (87-68-3)	
BCF fish 1	17000 (Salmo gairdneri)
BCF fish 2	7000 (Pleuronectes platessa, Flow-through system)
BCF other aquatic organisms 1	45.36 (Procambarus sp., Flow-through system)
BCF other aquatic organisms 2	3000 (Mytilus edulis, Flow-through system)
Log Pow	3.74 - 4.90
Bioaccumulative potential	High potential for bioaccumulation (BCF > 5000).
hexachloroethane (67-72-1)	
BCF fish 1	1200 (Salmo gairdneri)
BCF fish 2	756 mg/l (768 h, Pimephales promelas)
Log Pow	3.34 - 4.62
Bioaccumulative potential	Potential for bioaccumulation ($500 \leq \text{BCF} \leq 5000$).
indeno(1,2,3-cd)pyrene (193-39-5)	
BCF other aquatic organisms 1	10000 (240 h, Amphipoda)
Log Pow	6.6 - 7.7
Bioaccumulative potential	High potential for bioaccumulation (BCF > 5000).
naphthalene (91-20-3)	
BCF fish 1	23 - 168 (BCF; 8 weeks; Cyprinus carpio)
Log Pow	3.3 (Experimental value)
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).
nitrobenzene (98-95-3)	
BCF fish 1	15 (BCF; 672 h)
BCF fish 2	1.6 - 7.7 (BCF; OECD 305: Bioconcentration: Flow-Through Fish Test; 42 days; Cyprinus carpio; Flow-through system; Fresh water; Experimental value)
BCF other aquatic organisms 1	24 (BCF)
Log Pow	1.85 (Calculated; 1.86; Experimental value; EU Method A.8: Partition Coefficient)
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).
N-Nitrosodimethylamine (62-75-9)	
Log Pow	-0.77 - -0.57
Bioaccumulative potential	Bioaccumulation: not applicable.
2,3,4,5,6-pentachlorophenol (87-86-5)	
BCF fish 1	770 (768 h, Pimephales promelas)
BCF fish 2	39 - 224 (Cyprinus carpio, Test duration: 8 weeks)
BCF other aquatic organisms 1	1250 (Algae)
Log Pow	4.07 - 5.19
Bioaccumulative potential	Potential for bioaccumulation ($500 \leq \text{BCF} \leq 5000$).
Methylene Chloride (75-09-2)	
BCF fish 1	2 - 40 (OECD 305: Bioconcentration: Flow-Through Fish Test, 6 week(s), Cyprinus carpio, Semi-static system, Fresh water, Experimental value, GLP)
Log Pow	1.25 (Experimental value, OECD 107: Partition Coefficient (n-octanol/water): Shake Flask Method, 20 °C)
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).

12.4. Mobility in soil

atrazine (1912-24-9)	
Ecology - soil	Toxic to flora. Not toxic to bees.
benzo[a]anthracene (56-55-3)	
Ecology - soil	Adsorbs into the soil.

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Benzo(b)fluoranthene (205-99-2)	
Ecology - soil	Adsorbs into the soil.
benzo[k]fluoranthene (207-08-9)	
Ecology - soil	Adsorbs into the soil.
benzo[a]pyrene (50-32-8)	
Ecology - soil	Adsorbs into the soil.
chrysene (218-01-9)	
Ecology - soil	Adsorbs into the soil.
dibenz(a,h)anthracene (53-70-3)	
Ecology - soil	Adsorbs into the soil.
Bis(2-ethylhexyl) phthalate (117-81-7)	
Surface tension	0.032 N/m (20 °C)
Log Koc	5.2 (log Koc, Calculated value)
Ecology - soil	Adsorbs into the soil. Low potential for mobility in soil.
hexachlorobenzene (118-74-1)	
Ecology - soil	Adsorbs into the soil. Not toxic to bees.
hexachlorobuta-1,3-diene (87-68-3)	
Ecology - soil	Soil contaminant.
indeno(1,2,3-cd)pyrene (193-39-5)	
Ecology - soil	Adsorbs into the soil.
naphthalene (91-20-3)	
Surface tension	0.03 N/m (100 °C)
nitrobenzene (98-95-3)	
Surface tension	0.0439 N/m
Log Koc	Koc,Other; 118; Calculated value; log Koc; Other; 2.07; Calculated value
Ecology - soil	Low potential for adsorption in soil.
2,3,4,5,6-pentachlorophenol (87-86-5)	
Ecology - soil	No (test)data on mobility of the substance available.
Methylene Chloride (75-09-2)	
Surface tension	0.028 N/m (20 °C)
Ecology - soil	Low potential for adsorption in soil. May be harmful to plant growth, blooming and fruit formation.

12.5. Other adverse effects

Custom 8270 Short List	
atrazine (1912-24-9)	
benzo[a]anthracene (56-55-3)	
Benzo(b)fluoranthene (205-99-2)	
benzo[k]fluoranthene (207-08-9)	
benzo[a]pyrene (50-32-8)	
chrysene (218-01-9)	
dibenz(a,h)anthracene (53-70-3)	

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Bis(2-ethylhexyl) phthalate (117-81-7)	
hexachlorobenzene (118-74-1)	
hexachlorobuta-1,3-diene (87-68-3)	
hexachloroethane (67-72-1)	
indeno(1,2,3-cd)pyrene (193-39-5)	
naphthalene (91-20-3)	
nitrobenzene (98-95-3)	
N-Nitrosodimethylamine (62-75-9)	
2,3,4,5,6-pentachlorophenol (87-86-5)	
Methylene Chloride (75-09-2)	

Other information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

Department of Transportation (DOT)

In accordance with DOT

Transport document description : UN1593 Dichloromethane, 6.1, III
UN-No.(DOT) : UN1593
Proper Shipping Name (DOT) : Dichloromethane
Class (DOT) : 6.1 - Class 6.1 - Poisonous materials 49 CFR 173.132
Packing group (DOT) : III - Minor Danger
Hazard labels (DOT) : 6.1 - Poison



DOT Packaging Non Bulk (49 CFR 173.xxx) : 203
DOT Packaging Bulk (49 CFR 173.xxx) : 241

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

DOT Special Provisions (49 CFR 172.102)	: IB3 - Authorized IBCs: Metal (31A, 31B and 31N); Rigid plastics (31H1 and 31H2); Composite (31HZ1 and 31HA2, 31HB2, 31HN2, 31HD2 and 31HH2). Additional Requirement: Only liquids with a vapor pressure less than or equal to 110 kPa at 50 C (1.1 bar at 122 F), or 130 kPa at 55 C (1.3 bar at 131 F) are authorized, except for UN2672 (also see Special Provision IP8 in Table 2 for UN2672). IP8 - Ammonia solutions may be transported in rigid or composite plastic IBCs (31H1, 31H2 and 31HZ1) that have successfully passed, without leakage or permanent deformation, the hydrostatic test specified in 178.814 of this subchapter at a test pressure that is not less than 1.5 times the vapor pressure of the contents at 55 C (131 F). N36 - Aluminum or aluminum alloy construction materials are permitted only for halogenated hydrocarbons that will not react with aluminum. T7 - 4 178.274(d)(2) Normal..... 178.275(d)(3) TP2 - a. The maximum degree of filling must not exceed the degree of filling determined by the following: (image) Where: tr is the maximum mean bulk temperature during transport, tf is the temperature in degrees celsius of the liquid during filling, and a is the mean coefficient of cubical expansion of the liquid between the mean temperature of the liquid during filling (tf) and the maximum mean bulk temperature during transportation (tr) both in degrees celsius. b. For liquids transported under ambient conditions may be calculated using the formula: (image) Where: d15 and d50 are the densities (in units of mass per unit volume) of the liquid at 15 C (59 F) and 50 C (122 F), respectively.
DOT Packaging Exceptions (49 CFR 173.xxx)	: 153
DOT Quantity Limitations Passenger aircraft/rail (49 CFR 173.27)	: 60 L
DOT Quantity Limitations Cargo aircraft only (49 CFR 175.75)	: 220 L
DOT Vessel Stowage Location	: A - The material may be stowed "on deck" or "under deck" on a cargo vessel and on a passenger vessel.
Emergency Response Guide (ERG) Number	: 160
Other information	: No supplementary information available.

Transportation of Dangerous Goods

Not applicable

Transport by sea

Transport document description (IMDG)	: UN 1593 DICHLOROMETHANE, 6.1, III, MARINE POLLUTANT/ENVIRONMENTALLY HAZARDOUS
UN-No. (IMDG)	: 1593
Proper Shipping Name (IMDG)	: DICHLOROMETHANE
Class (IMDG)	: 6.1 - Toxic substances
Packing group (IMDG)	: III - substances presenting low danger

Air transport

Transport document description (IATA)	: UN 1593 Dichloromethane, 6.1, III, ENVIRONMENTALLY HAZARDOUS
UN-No. (IATA)	: 1593
Proper Shipping Name (IATA)	: Dichloromethane
Class (IATA)	: 6.1 - Toxic Substances
Packing group (IATA)	: III - Minor Danger

SECTION 15: Regulatory information

15.1. US Federal regulations

atrazine (1912-24-9)

Listed on the United States TSCA (Toxic Substances Control Act) inventory
Subject to reporting requirements of United States SARA Section 313

benzo[a]anthracene (56-55-3)

Listed on the United States TSCA (Toxic Substances Control Act) inventory
Subject to reporting requirements of United States SARA Section 313

CERCLA RQ	10 lb
-----------	-------

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Benzo(b)fluoranthene (205-99-2)	
Not listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	1 lb
benzo[k]fluoranthene (207-08-9)	
Not listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	5000 lb
benzo[a]pyrene (50-32-8)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	1 lb
chrysene (218-01-9)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	100 lb
dibenz(a,h)anthracene (53-70-3)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	1 lb
Bis(2-ethylhexyl) phthalate (117-81-7)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	100 lb
hexachlorobenzene (118-74-1)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	10 lb
hexachlorobuta-1,3-diene (87-68-3)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	1 lb
hexachloroethane (67-72-1)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
EPA TSCA Regulatory Flag	TP - TP - indicates a substance that is the subject of a proposed TSCA section 4 test rule.
CERCLA RQ	100 lb
indeno(1,2,3-cd)pyrene (193-39-5)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
CERCLA RQ	100 lb
naphthalene (91-20-3)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	100 lb

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

nitrobenzene (98-95-3)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	1000 lb
RQ (Reportable quantity, section 304 of EPA's List of Lists)	1000 lb
SARA Section 302 Threshold Planning Quantity (TPQ)	10000 lb
N-Nitrosodimethylamine (62-75-9)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	10 lb
RQ (Reportable quantity, section 304 of EPA's List of Lists)	10 lb
SARA Section 302 Threshold Planning Quantity (TPQ)	1000 lb
2,3,4,5,6-pentachlorophenol (87-86-5)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
CERCLA RQ	10 lb
Methylene Chloride (75-09-2)	
Listed on the United States TSCA (Toxic Substances Control Act) inventory Subject to reporting requirements of United States SARA Section 313	
Listed on EPA Hazardous Air Pollutant (HAPS)	
EPA TSCA Regulatory Flag	R - R - indicates a substance that is the subject of a TSCA section 6 risk management rule.
CERCLA RQ	1000 lb

15.2. International regulations

CANADA

atrazine (1912-24-9)	
Listed on the Canadian DSL (Domestic Substances List)	
benzo[a]anthracene (56-55-3)	
Listed on the Canadian NDSL (Non-Domestic Substances List)	
Benzo(b)fluoranthene (205-99-2)	
Not listed on the Canadian DSL (Domestic Substances List)/NDSL (Non-Domestic Substances List)	
benzo[k]fluoranthene (207-08-9)	
Not listed on the Canadian DSL (Domestic Substances List)/NDSL (Non-Domestic Substances List)	
benzo[a]pyrene (50-32-8)	
Listed on the Canadian DSL (Domestic Substances List)	
chrysene (218-01-9)	
Listed on the Canadian DSL (Domestic Substances List)	
dibenz(a,h)anthracene (53-70-3)	
Listed on the Canadian NDSL (Non-Domestic Substances List)	
Bis(2-ethylhexyl) phthalate (117-81-7)	
Listed on the Canadian DSL (Domestic Substances List)	
hexachlorobenzene (118-74-1)	
Listed on the Canadian DSL (Domestic Substances List)	
hexachlorobuta-1,3-diene (87-68-3)	
Listed on the Canadian DSL (Domestic Substances List)	
hexachloroethane (67-72-1)	
Listed on the Canadian DSL (Domestic Substances List)	

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

indeno(1,2,3-cd)pyrene (193-39-5)

Listed on the Canadian NDSL (Non-Domestic Substances List)

naphthalene (91-20-3)

Listed on the Canadian DSL (Domestic Substances List)

nitrobenzene (98-95-3)

Listed on the Canadian DSL (Domestic Substances List)

N-Nitrosodimethylamine (62-75-9)

Listed on the Canadian NDSL (Non-Domestic Substances List)

2,3,4,5,6-pentachlorophenol (87-86-5)

Listed on the Canadian DSL (Domestic Substances List)

Methylene Chloride (75-09-2)

Listed on the Canadian DSL (Domestic Substances List)

EU-Regulations

No additional information available

National regulations

benzo[a]anthracene (56-55-3)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

Benzo(b)fluoranthene (205-99-2)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

benzo[k]fluoranthene (207-08-9)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

benzo[a]pyrene (50-32-8)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

chrysene (218-01-9)

Listed on IARC (International Agency for Research on Cancer)

dibenz(a,h)anthracene (53-70-3)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

Bis(2-ethylhexyl) phthalate (117-81-7)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)
Listed on EPA Hazardous Air Pollutant (HAPS)

hexachlorobenzene (118-74-1)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)
Listed on EPA Hazardous Air Pollutant (HAPS)

hexachlorobuta-1,3-diene (87-68-3)

Listed on EPA Hazardous Air Pollutant (HAPS)

hexachloroethane (67-72-1)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)
Listed on EPA Hazardous Air Pollutant (HAPS)

indeno(1,2,3-cd)pyrene (193-39-5)

Listed on IARC (International Agency for Research on Cancer)
Listed as carcinogen on NTP (National Toxicology Program)

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

naphthalene (91-20-3)
Listed on IARC (International Agency for Research on Cancer) Listed as carcinogen on NTP (National Toxicology Program) Listed on EPA Hazardous Air Pollutant (HAPS)
nitrobenzene (98-95-3)
Listed on IARC (International Agency for Research on Cancer) Listed as carcinogen on NTP (National Toxicology Program) Listed on EPA Hazardous Air Pollutant (HAPS)
N-Nitrosodimethylamine (62-75-9)
Listed on IARC (International Agency for Research on Cancer) Listed as carcinogen on NTP (National Toxicology Program) Listed on EPA Hazardous Air Pollutant (HAPS)
2,3,4,5,6-pentachlorophenol (87-86-5)
Listed on IARC (International Agency for Research on Cancer) Listed as carcinogen on NTP (National Toxicology Program) Listed on EPA Hazardous Air Pollutant (HAPS)
Methylene Chloride (75-09-2)
Listed on IARC (International Agency for Research on Cancer) Listed as carcinogen on NTP (National Toxicology Program) Listed on EPA Hazardous Air Pollutant (HAPS)

15.3. US State regulations

atrazine (1912-24-9)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
No	Yes	Yes	Yes		100 µg/day (oral)
benzo[a]anthracene (56-55-3)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.033 µg/day	
Benzo(b)fluoranthene (205-99-2)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.096 µg/day	
benzo[k]fluoranthene (207-08-9)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No		
benzo[a]pyrene (50-32-8)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.06 µg/day	

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

chrysene (218-01-9)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.35 µg/day	
dibenz(a,h)anthracene (53-70-3)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.2 µg/day	
Bis(2-ethylhexyl) phthalate (117-81-7)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	Yes	No	Yes	310 µg/day	4200 µg/day (intravenous), Adult; 600 µg/day (intravenous), Infant boys, age 29 days - 24 mos; 210 µg/day (intravenous), Neonatal infant boys, age 0 - 28 days; 410 µg/day (oral), Adult; 58 µg/day (oral), Infant boys, age 29 days - 24 mos; 20 µg/day (oral), Neonatal infant boys, age 0 - 28 days
hexachlorobenzene (118-74-1)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	Yes	No	No	0.4 µg/day	
hexachlorobuta-1,3-diene (87-68-3)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No		
hexachloroethane (67-72-1)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	20 µg/day	

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

indeno(1,2,3-cd)pyrene (193-39-5)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No		
naphthalene (91-20-3)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	5.8 µg/day	
nitrobenzene (98-95-3)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	Yes		
N-Nitrosodimethylamine (62-75-9)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	0.04 µg/day	
2,3,4,5,6-pentachlorophenol (87-86-5)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	40 µg/day	
Methylene Chloride (75-09-2)					
U.S. - California - Proposition 65 - Carcinogens List	U.S. - California - Proposition 65 - Developmental Toxicity	U.S. - California - Proposition 65 - Reproductive Toxicity - Female	U.S. - California - Proposition 65 - Reproductive Toxicity - Male	No significant risk level (NSRL)	Maximum allowable dose level (MADL)
Yes	No	No	No	50 µg/day	

SECTION 16: Other information

Revision date : 08/07/2019

Other information : None.

Full text of H-phrases:

H227	Combustible liquid
H317	May cause an allergic skin reaction
H340	May cause genetic defects
H350	May cause cancer
H360	May damage fertility or the unborn child

Custom 8270 Short List

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.