

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Product form : Mixture  
 Product name : roQ™ QuEChERS dSPE Kit - 2mL CT, 150mg MgSO<sub>4</sub>, 25mg PSA, 25mg C18E  
 Product code : KS0-8913

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

### 1.3. Details of the supplier of the safety data sheet

Phenomenex, Inc.  
 411 Madrid Ave.  
 90501 Torrance - United States  
 T 310-212-0555 - F 310-328-7768  
[info@phenomenex.com](mailto:info@phenomenex.com) - [www.phenomenex.com](http://www.phenomenex.com)

### 1.4. Emergency telephone number

Emergency number : INFORTAC assistance (United States): +1-800-535-5053  
 INFOTRAC assistance (Outside United States): +1-352-323-3500

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Classification (GHS-US)

Not classified

### 2.2. Label elements

#### GHS-US labelling

No labelling applicable

### 2.3. Other hazards

No additional information available

### 2.4. Unknown acute toxicity (GHS US)

No data available

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Not applicable

### 3.2. Mixtures

Name	Product identifier	%	Classification (GHS-US)
magnesium sulfate (Component)	(CAS No.) 7487-88-9	75	
Sepra C18E Bonded Silica (Component)	(CAS No.) 112926-00-8	12.5	
Sepra PSA Bonded Silica (Component)	(CAS No.) 112926-00-8	12.5	

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general : Consult a physician after significant exposure.  
 First-aid measures after inhalation : Remove victim to fresh air. If not breathing, give artificial respiration, preferably mouth to mouth.  
 First-aid measures after skin contact : Wash immediately with plenty of soap and water.  
 First-aid measures after eye contact : Rinse immediately with plenty of water.  
 First-aid measures after ingestion : Rinse mouth out with water. Never give anything by mouth to an unconscious person.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire.

### 5.2. Special hazards arising from the substance or mixture

Fire hazard : No fire hazard.

# roQ™ QuEChERS dSPE Kit - 2mL CT, 150mg MgSO<sub>4</sub>, 25mg PSA, 25mg C18E

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Explosion hazard : No direct explosion hazard.  
Reactivity : No data available.

### 5.3. Advice for firefighters

Firefighting instructions : No specific measures are necessary.  
Protection during firefighting : Wear a self contained breathing apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Protective equipment : Use personal protective equipment as required. Avoid breathing dust, mist or spray.

#### 6.1.2. For emergency responders

Protective equipment : Use personal protective equipment as required. Avoid breathing dust, mist or spray.

### 6.2. Environmental precautions

Do not allow into drains or water courses or dispose of where ground or surface waters may be affected.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Large spills: scoop solid spill into closing containers.

### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Use personal protective equipment as required. Local exhaust is recommended where dust may occur.

Hygiene measures : Handle in accordance with good industrial hygiene and safety procedures.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container tightly closed and in well ventilated place.

### 7.3. Specific end use(s)

No additional information available


## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

### 8.2. Exposure controls

Appropriate engineering controls : Provide local exhaust or general room ventilation to minimize exposure to dust.

Personal protective equipment : Gloves. Safety glasses.



Hand protection : Wear suitable gloves resistant to chemical penetration.

Eye protection : Wear closed safety glasses.

Skin and body protection : Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin contact.

Respiratory protection : Respiratory protection not required in normal conditions.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Solid  
Colour : No data available  
Odour : No data available  
Odour threshold : No data available  
pH : No data available  
Relative evaporation rate (butylacetate=1) : No data available  
Melting point : No data available  
Freezing point : No data available  
Boiling point : No data available  
Flash point : No data available

# roQ™ QuEChERS dSPE Kit - 2mL CT, 150mg MgSO4, 25mg PSA, 25mg C18E

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Self ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Log Pow	: No data available
Log Kow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No data available.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No additional information available

### 10.4. Conditions to avoid

No additional information available

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
Specific target organ toxicity (single exposure)	: Not classified
Specific target organ toxicity (repeated exposure)	: Not classified
Aspiration hazard	: Not classified

## SECTION 12: Ecological information

### 12.1. Toxicity

No additional information available

### 12.2. Persistence and degradability

No additional information available

### 12.3. Bioaccumulative potential

No additional information available

### 12.4. Mobility in soil

No additional information available

# roQ™ QuEChERS dSPE Kit - 2mL CT, 150mg MgSO4, 25mg PSA, 25mg C18E

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### 12.5. Other adverse effects

No additional information available

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Waste disposal recommendations : Consult an expert on waste disposal or treatment. Dispose of at a licensed waste collection center.

## SECTION 14: Transport information

In accordance with ADR / RID / ADN / IMDG / ICAO / IATA

### 14.1. UN number

Not applicable

### 14.2. UN proper shipping name

Not applicable

### 14.3. Additional information

#### Overland transport

No additional information available

#### Transport by sea

No additional information available

#### Air transport

No additional information available

## SECTION 15: Regulatory information

### 15.1. US Federal regulations

No additional information available

### 15.2. International regulations

#### CANADA

No additional information available

#### EU-Regulations

No additional information available

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

#### Classification according to Directive 67/548/EEC or 1999/45/EC

Not classified

### 15.2.2. National regulations

No additional information available

### 15.3. US State regulations

No additional information available


# roQ™ QuEChERS dSPE Kit - 2mL CT, 150mg MgSO4, 25mg PSA, 25mg C18E

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### SECTION 16: Other information

- NFPA health hazard : 0 - Exposure under fire conditions would offer no hazard beyond that of ordinary combustible materials.
- NFPA fire hazard : 0 - Materials that will not burn.
- NFPA reactivity : 1 - Normally stable, but can become unstable at elevated temperatures and pressures or may react with water with some release of energy, but not violently.



#### HMIS III Rating

- Health : 0 Minimal Hazard - No significant risk to health
- Flammability : 0 Minimal Hazard
- Physical : 0 Minimal Hazard
- Personal Protection : B

PHX US GHS

Copyright 2012 Phenomenex, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the the properties of the product. Phenomenex, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at [www.phenomenex.com](http://www.phenomenex.com) for additional terms and conditions of sale.