

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

Date of issue: 22/10/2015 Revision date: : Version: 1.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : o-Terphenyl Standard

Product code : AL0-130012
Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Laboratory Use Industrial/Professional use spec : Industrial

For professional use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Phenova

6390 Joyce Dr. Suite 100

80403 Golden, CO - United States T 1-866-942-2978 - F 1-866-283-0269

info@phenova.com - www.phenova.com

1.4. Emergency telephone number

Emergency number : ChemTel Assistance (US/Canada) 1-800-255-3924

ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Flam. Liq. 2 H225 STOT SE 3 H336 Aquatic Chronic 2 H411

Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]

F+; R12 N; R51/53 R66 R67

Full text of R-phrases: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :







GHS07

GHS09

Signal word (CLP) : Danger Hazardous ingredients : n-pentane

Hazard statements (CLP) : H225 - Highly flammable liquid and vapor

H336 - May cause drowsiness or dizziness H411 - Toxic to aquatic life with long lasting effects

27/10/2015 EN (English US) 1/7

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

Precautionary statements (CLP) : P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking

P233 - Keep container tightly closed

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray P271 - Use only outdoors or in a well-ventilated area

P273 - Avoid release to the environment

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P308+P313 - IF exposed or concerned: Get medical advice/attention

P391 - Collect spillage

P403+P233 - Store in a well-ventilated place. Keep container tightly closed

EUH phrases : EUH066 - Repeated exposure may cause skin dryness or cracking

No labeling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
n-pentane	(CAS No) 109-66-0 (EC no) 203-692-4 (EC index no) 601-006-00-1	95	Flam. Liq. 2, H225 STOT SE 3, H336 Asp. Tox. 1, H304 Aquatic Chronic 2, H411
o-terphenyl (Component)	(CAS No) 84-15-1 (EC no) 201-517-6	5	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319 STOT SE 3, H335

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. IF exposed or concerned: Get

medical advice/attention.

First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.

First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed

by warm water rinse

First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness

persist.

First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

No additional information available

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

27/10/2015 EN (English US) 2/7

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.

Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect

spillage. Store away from other materials.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or

smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. Obtain special instructions before use. Do not handle until all safety precautions have

been read and understood.

Hygiene measures : Gently wash with plenty of soap and water. Remove/Take off immediately all contaminated

clothing. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container closed when not in use. Keep container tightly closed and in a well-ventilated

place. Keep away from any flames or sparking source.

Incompatible products : Strong bases. Strong acids.

Incompatible materials : Sources of ignition. Direct sunlight.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Appropriate engineering controls : Either local exhaust or general room ventilation is usually required.

Personal protective equipment : Avoid all unnecessary exposure. Gloves. Protective clothing. Protective goggles. Safety

glasses







Hand protection : Wear chemically resistant protective gloves. Wear suitable gloves resistant to chemical

penetration.

Eye protection : Chemical goggles or safety glasses. Safety glasses.

Skin and body protection : Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin

contact.

Respiratory protection : Wear appropriate mask.

Other information : Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Liquid Color Colorless Odor characteristic рΗ No data available Melting point No data available Freezing point No data available Boiling point No data available No data available Flash point Auto-ignition temperature No data available Decomposition temperature : No data available

27/10/2015 EN (English US) 3/7

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

Flammability (solid, gas) : Non flammable.

Relative density : No data available

Solubility : No data available

Explosive properties : No data available

Oxidizing properties : No data available

Explosion limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Not established.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

fume. Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

o-terphenyl (84-15-1)		
LD50 oral rat	2000 mg/kg (Rat)	
ATE CLP (oral)	2000.000 mg/kg body weight	
n-pentane (109-66-0)		
LD50 oral rat	> 2000 mg/kg (Rat; OECD 401: Acute Oral Toxicity; Experimental value)	

Skin corrosion/irritation : Not classified

Based on available data, the classification criteria are not met

Serious eye damage/irritation : Not classified

Based on available data, the classification criteria are not met

Respiratory or skin sensitization : Not classified

Based on available data, the classification criteria are not met

Germ cell mutagenicity : Not classified

Based on available data, the classification criteria are not met

Carcinogenicity : Not classified

May cause cancer

Reproductive toxicity : Not classified

Based on available data, the classification criteria are not met

Specific target organ toxicity (single exposure) : May cause drowsiness or dizziness.

Based on available data, the classification criteria are not met

Specific target organ toxicity (repeated

exposure)

: Not classified

Based on available data, the classification criteria are not met

Aspiration hazard : Not classified

Based on available data, the classification criteria are not met : Based on available data, the classification criteria are not met.

Potential Adverse human health effects and

symptoms

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

27/10/2015 EN (English US) 4/7

Class (IATA) Class (IMDG)

Class (ADN)

o-Terphenyl Standard
Safety Data Sheet
according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

	1) with its amendment Regulation (EC) No. 453/2010	
12.2. Persistence and degradability		
o-Terphenyl Standard		
Persistence and degradability	Not established.	
o-terphenyl (84-15-1)		
Persistence and degradability	Not readily biodegradable in water. Forming sediments in water. Biodegradability in soil: no data available. Adsorbs into the soil.	
ThOD	2.99 g O /g substance	
n-pentane (109-66-0)		
Persistence and degradability	Readily biodegradable in water. Low potential for adsorption in soil.	
12.3. Bioaccumulative potential		
o-Terphenyl Standard		
Bioaccumulative potential	Not established.	
o-terphenyl (84-15-1)		
BCF fish 1	590 (BCF)	
Log Pow	5.52 (Estimated value)	
Bioaccumulative potential	High potential for bioaccumulation (Log Kow > 5).	
n-pentane (109-66-0)		
BCF fish 1	171 (BCF)	
Log Pow	3.45 (Experimental value; 25 °C)	
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).	
12.4. Mobility in soil		
n-pentane (109-66-0)		
Surface tension	0.015 N/m (25 °C; 100 %; 0.013 N/m; 20 °C)	
Log Koc	log Koc,2.9; QSAR	
12.5. Results of PBT and vPvB assess	ment	
No additional information available		
12.6. Other adverse effects		
Additional information	: Avoid release to the environment	
SECTION 13: Disposal considerat	ions	
13.1. Waste treatment methods		
Waste disposal recommendations	: Dispose in a safe manner in accordance with local/national regulations.	
Ecology - waste materials	: Avoid release to the environment.	
SECTION 14: Transport information	on	
In accordance with ADR / RID / IMDG / IATA		
14.1. UN number		
UN-No. (ADR)	: 1993	
UN-No.(IATA)	: 1993	
14.2. UN proper shipping name		
Proper Shipping Name (ADR)	: FLAMMABLE LIQUID, N.O.S.	
Proper Shipping Name (IATA)	: FLAMMABLE LIQUID, N.O.S.	
Proper Shipping Name (IMDG)	: FLAMMABLE LIQUID, N.O.S.	
Proper Shipping Name (ADN)	: FLAMMABLE LIQUID, N.O.S.	
Transport document description (ADR)	: UN 1993 FLAMMABLE LIQUID, N.O.S. (dichloromethane(75-09-2)), 3, II, (D/E),	
	ENVIRONMENTALLY HAZARDOUS	
14.3. Packing group		
Class (ADR)	: 3	
Classification code (ADR)	: F1	

27/10/2015 EN (English US) 5/7

: 3

: 6.1

: 6.1

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

Hazard labels (ADR) : 3



Hazard labels (IATA) : 3



14.4. Packing group

Packing group (ADR) : II Packing group (IATA) : II

14.5. Environmental hazards

Dangerous for the environment



Other information : No supplementary information available.

14.6. Special precautions for user

14.6.1. Overland transport

Hazard identification number (Kemler No.) : 33 Classification code (ADR) : F1

Orange plates :

33 1993

Special provision (ADR) : 274, 601, 640D

Transport category (ADR) : 2
Tunnel restriction code (ADR) : D/E
Limited quantities (ADR) : 11
Excepted quantities (ADR) : E2

14.6.2. Transport by sea

No additional information available

14.6.3. Air transport

CAO packing instructions (IATA) : 364 CAO max net quantity (IATA) : 60L PCA packing instructions (IATA) : 353 PCA Limited quantities (IATA) : Y341 PCA limited quantity max net quantity (IATA) : 1L : 5L PCA max net quantity (IATA) PCA Excepted quantities (IATA) : E2 ERG code (IATA) : 3H

14.6.4. Inland waterway transport

Carriage prohibited (ADN) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

27/10/2015 EN (English US) 6/7

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EC) No. 453/2010

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending

Regulation (EC) No 1907/2006.

Other information : None.

PHV SDS EU

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.

27/10/2015 EN (English US) 7/7