

Safety Data Sheet

Date of issue: 29/02/2016 Revision date: : Version: 1.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : 8141 OP Pesticides Mix 2

Product code : AL0-101541
Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Laboratory Use Industrial/Professional use spec : Industrial

For professional use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Phenova

6390 Joyce Dr. Suite 100

Emergency number

80403 Golden, CO - United States T 1-866-942-2978 - F 1-866-283-0269

info@phenova.com - www.phenova.com

1.4. Emergency telephone number

: ChemTel Assistance (US/Canada) 1-800-255-3924 ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Flam. Liq. 2 H225
Skin Irrit. 2 H315
Eye Irrit. 2 H319
Repr. 2 H361
STOT SE 3 H336
STOT RE 2 H373
Aquatic Acute 1 H400
Aquatic Chronic 2 H411

Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]

Repr.Cat.3; R62 F; R11 Xn; R48/20 Xi; R38 N; R51/53 R67

Full text of R-phrases: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

29/02/2016 EN (English US) 1/9

Safety Data Sheet

2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)









Signal word (CLP) : Danger Hazardous ingredients : hexane

Hazard statements (CLP) : H225 - Highly flammable liquid and vapor

H315 - Causes skin irritation

H319 - Causes serious eye irritation

H336 - May cause drowsiness or dizziness

H361 - Suspected of damaging fertility or the unborn child

H373 - May cause damage to organs through prolonged or repeated exposure

H410 - Very toxic to aquatic life with long lasting effects

Precautionary statements (CLP) : P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking

P233 - Keep container tightly closed

P260 - Do not breathe dust/fume/gas/mist/vapors/spray P271 - Use only outdoors or in a well-ventilated area

P273 - Avoid release to the environment

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P302+P350 - IF ON SKIN: Gently wash with plenty of soap and water

P304+P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing

P308+P313 - IF exposed or concerned: Get medical advice/attention P337+P313 - If eye irritation persists: Get medical advice/attention P362+P364 - Take off contaminated clothing and wash it before reuse

P391 - Collect spillage

P403+P235 - Store in a well-ventilated place. Keep cool

No labeling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
hexane (Component)	(CAS No) 110-54-3 (EC no) 203-777-6 (EC index no) 601-037-00-0	89.86	Flam. Liq. 2, H225 Skin Irrit. 2, H315 Repr. 2, H361f STOT SE 3, H336 STOT RE 2, H373 Asp. Tox. 1, H304 Aquatic Chronic 2, H411
acetone (Component)	(CAS No) 67-64-1 (EC no) 200-662-2 (EC index no) 606-001-00-8	10	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336
sulfotep (Component)	(CAS No) 3689-24-5 (EC no) 222-995-2 (EC index no) 015-027-00-3	0.02	Acute Tox. 1 (Oral), H300 Acute Tox. 1 (Dermal), H310 Aquatic Acute 1, H400 (M=1000) Aquatic Chronic 1, H410
parathion (Component)	(CAS No) 56-38-2 (EC no) 200-271-7 (EC index no) 015-034-00-1	0.02	Acute Tox. 1 (Oral), H300 Acute Tox. 1 (Dermal), H310 Acute Tox. 2 (Inhalation), H330 STOT RE 1, H372 Aquatic Acute 1, H400 (M=100) Aquatic Chronic 1, H410
malathion (Component)	(CAS No) 121-75-5 (EC no) 204-497-7 (EC index no) 015-041-00-X	0.02	Acute Tox. 3 (Oral), H301 Skin Sens. 1, H317 Aquatic Acute 1, H400 (M=1000) Aquatic Chronic 1, H410

29/02/2016 EN (English US) 2/9

Safety Data Sheet

Name	Product identifier	Specific concentration limits
hexane (Component)	(CAS No) 110-54-3 (EC no) 203-777-6 (EC index no) 601-037-00-0	(C >= 5) STOT RE 2, H373

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. IF exposed or concerned: Get

medical advice/attention.

First-aid measures after inhalation : Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a

POISON CENTER or doctor/physician if you feel unwell.

First-aid measures after skin contact : Rinse skin with water/shower. Remove/Take off immediately all contaminated clothing. Wash

with plenty of soap and water. Wash contaminated clothing before reuse. If skin irritation

occurs: Get medical advice/attention.

First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness

persist.

First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries after inhalation : May cause drowsiness or dizziness.

Symptoms/injuries after skin contact : Causes skin irritation.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Highly flammable liquid and vapor.

Explosion hazard : May form flammable/explosive vapor-air mixture.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection. Avoid breathing dust/fume/gas/mist/vapors/spray.

Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up in absorbent material. Collect spillage.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Additional hazards when processed

: Handle empty containers with care because residual vapors are flammable.

Precautions for safe handling

: Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. No open flames. No smoking. Use only non-sparking tools. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Use only

outdoors or in a well-ventilated area.

Hygiene measures : Gently wash with plenty of soap and water. Remove/Take off immediately all contaminated

clothing. Wash contaminated clothing before reuse.

29/02/2016 EN (English US) 3/9

Safety Data Sheet

7.2. Conditions for safe storage, including any incompatibilities

Technical measures : Proper grounding procedures to avoid static electricity should be followed. Ground/bond

container and receiving equipment.

Storage conditions : Keep in fireproof place. Keep container tightly closed. Keep container tightly closed and in a

well-ventilated place. Keep away from any flames or sparking source.

Incompatible materials : Direct sunlight. Heat sources.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Appropriate engineering controls : Either local exhaust or general room ventilation is usually required.

Personal protective equipment : Avoid all unnecessary exposure. Gloves. Protective clothing. Protective goggles. Safety

glasses.







Hand protection : Wear chemically resistant protective gloves. Wear suitable gloves resistant to chemical

penetration.

Eye protection : Chemical goggles or safety glasses. Safety glasses.

Skin and body protection : Wear suitable protective clothing. Wear chemically protective gloves, lab coat or apron to

prevent prolonged or repeated skin contact.

Respiratory protection : Where exposure through inhalation may occur from use, respiratory protection equipment is

recommended.

Other information : Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Liquid Color Colorless. Odor characteristic. рΗ No data available Melting point No data available Freezing point No data available No data available Boiling point Flash point No data available Auto-ignition temperature No data available Decomposition temperature No data available

Flammability (solid, gas) : Highly flammable liquid and vapor

Relative density : No data available
Solubility : No data available
Explosive properties : No data available
Oxidizing properties : No data available
Explosion limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Highly flammable liquid and vapor. May form flammable/explosive vapor-air mixture.

10.3. Possibility of hazardous reactions

Not established.

29/02/2016 EN (English US) 4/9

Safety Data Sheet

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures. Open flame.

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

May release flammable gases.

SECTION 11: Toxicological information

11 1	Information on	

11.1. Information on toxicological effects	
Acute toxicity	: Not classified
sulfotep (3689-24-5)	
LD50 oral rat	5 mg/kg (Rat)
LD50 dermal rabbit	20 mg/kg (Rabbit)
ATE CLP (oral)	5.000 mg/kg body weight
ATE CLP (dermal)	20.000 mg/kg body weight
malathion (121-75-5)	
LD50 oral rat	290 mg/kg (Rat)
LD50 dermal rat	4444 mg/kg (Rat)
LD50 dermal rabbit	4100 mg/kg (Rabbit)
ATE CLP (oral)	290.000 mg/kg body weight
ATE CLP (dermal)	4100.000 mg/kg body weight
parathion (56-38-2)	
LD50 oral rat	2 mg/kg (Rat)
LD50 dermal rat	73 mg/kg (Rat)
LD50 dermal rabbit	40 mg/kg (Rabbit)
LC50 inhalation rat (mg/l)	0.03 mg/l/4h (Rat)
ATE CLP (oral)	2.000 mg/kg body weight
ATE CLP (dermal)	40.000 mg/kg body weight
ATE CLP (gases)	100.000 ppmV/4h
ATE CLP (vapors)	0.030 mg/l/4h
ATE CLP (dust, mist)	0.030 mg/l/4h
acetone (67-64-1)	
LD50 oral rat	5800 mg/kg (Rat; Equivalent or similar to OECD 401; Experimental value)
LD50 dermal rabbit	20000 mg/kg (Rabbit; Experimental value; Equivalent or similar to OECD 402; >7426 mg/kg bodyweight; Rabbit; Weight of evidence)
LC50 inhalation rat (mg/l)	71 mg/l/4h (Rat; Experimental value; 76 mg/l/4h; Rat; Experimental value)
LC50 inhalation rat (ppm)	30000 ppm/4h (Rat; Experimental value)
ATE CLP (oral)	5800.000 mg/kg body weight
ATE CLP (dermal)	20000.000 mg/kg body weight
ATE CLP (gases)	30000.000 ppmV/4h
ATE CLP (vapors)	71.000 mg/l/4h
ATE CLP (dust, mist)	71.000 mg/l/4h

hexane (110-54-3)	
LD50 oral rat	16000 mg/kg body weight (Rat; Equivalent or similar to OECD 401; Experimental value)
LD50 dermal rabbit	> 3350 mg/kg body weight (Rabbit; Read-across; Equivalent or similar to OECD 402)
ATE CLP (oral)	16000.000 mg/kg body weight

 ${\bf Skin\ corrosion/irritation}\qquad \qquad :\ {\bf Causes\ skin\ irritation}.$

Serious eye damage/irritation : Causes serious eye irritation.

Based on available data, the classification criteria are not met

Respiratory or skin sensitization : Not classified

Based on available data, the classification criteria are not met

Germ cell mutagenicity : Not classified

Based on available data, the classification criteria are not met

Carcinogenicity : Not classified

Based on available data, the classification criteria are not met

May cause cancer

Reproductive toxicity : Suspected of damaging fertility or the unborn child.

Specific target organ toxicity (single exposure) : May cause drowsiness or dizziness.

29/02/2016 EN (English US) 5/9

Safety Data Sheet

Specific target organ toxicity (repeated exposure)

: May cause damage to organs through prolonged or repeated exposure.

: Not classified Aspiration hazard

Based on available data, the classification criteria are not met

Potential Adverse human health effects and

symptoms

: Based on available data, the classification criteria are not met.

SECTION 12: Ecological in	normation
2.1. Toxicity	
cology - water	: Very toxic to aquatic life. Toxic to aquatic life with long lasting effects.
sulfotep (3689-24-5)	
EC50 Daphnia 1	0.00023 mg/l (EC50; 48 h)
LC50 fish 2	0.0036 mg/l (LC50; 96 h)
malathion (121-75-5)	
EC50 Daphnia 1	0.0008 mg/l (EC50; 48 h)
LC50 fish 2	0.17 mg/l (LC50; 96 h)
parathion (56-38-2)	
EC50 Daphnia 1	0.0025 mg/l (EC50; OECD 202: Daphnia sp. Acute Immobilisation Test; 48 h; Daphnia magna
LC50 fish 2	0.75 mg/l (LC50; 96 h)
acetone (67-64-1)	
LC50 fish 2	5540 mg/l (LC50; EU Method C.1; 96 h; Salmo gairdneri; Static system; Fresh water; Experimental value)
EC50 Daphnia 2	12600 mg/l (LC50; Other; 48 h; Daphnia magna; Static system; Fresh water; Experimental value)
hexane (110-54-3)	
LC50 fish 1	2.5 mg/l (LC50; 96 h)
EC50 Daphnia 1	2.1 mg/l (EC50; 48 h)
Threshold limit algae 2	26 mg/l (EbC50; OECD 201: Alga, Growth Inhibition Test; 72 h; Pseudokirchneriella subcapitata; Static system)

12.2. Persistence and degradability	
8141 OP Pesticides Mix 2	
Persistence and degradability	May cause long-term adverse effects in the environment.
malathion (121-75-5)	
Persistence and degradability	Biodegradable in the soil.
parathion (56-38-2)	
Persistence and degradability	Biodegradable in water. Forming sediments in water. Biodegradable in the soil. Adsorbs into the soil.
acetone (67-64-1)	
Persistence and degradability	Readily biodegradable in water. Biodegradable in the soil. Biodegradable in the soil under anaerobic conditions. No (test)data on mobility of the substance available.
Biochemical oxygen demand (BOD)	1.43 g O□ /g substance
Chemical oxygen demand (COD)	1.92 g O□ /g substance
ThOD	2.20 g O□ /g substance
BOD (% of ThOD)	0.872 (20 days; Literature study)
hexane (110-54-3)	
Persistence and degradability	Readily biodegradable in water. Photooxidation in water. easily degradable in the soil.
ThOD	3.52 g O□ /g substance
BOD (% of ThOD)	0.63 (Literature study)
12.3. Bioaccumulative potential	
8141 OP Pesticides Mix 2	
Bioaccumulative potential	Not established.
sulfotep (3689-24-5)	
Log Pow	3.99 (Experimental value)
Bioaccumulative potential	Bioaccumable.
malathion (121-75-5)	
Log Pow	2.36 - 2.89
Bioaccumulative potential	Low potential for bioaccumulation (Log Kow < 4).

29/02/2016 EN (English US) 6/9

335 (BCF; 912 h) 462 (BCF; 72 h) 240 (BCF; 999 h) 97 (BCF; 792 h) 3.8 Low potential for bioaccumulation (BCF < 500).
462 (BCF; 72 h) 240 (BCF; 999 h) 97 (BCF; 792 h) 3.8 Low potential for bioaccumulation (BCF < 500).
240 (BCF; 999 h) 97 (BCF; 792 h) 3.8 Low potential for bioaccumulation (BCF < 500). 0.69 (BCF) 3 (BCF; BCFWIN) -0.24 (Test data)
97 (BCF; 792 h) 3.8 Low potential for bioaccumulation (BCF < 500). 0.69 (BCF) 3 (BCF; BCFWIN) -0.24 (Test data)
3.8 Low potential for bioaccumulation (BCF < 500). 0.69 (BCF) 3 (BCF; BCFWIN) -0.24 (Test data)
Low potential for bioaccumulation (BCF < 500). 0.69 (BCF) 3 (BCF; BCFWIN) -0.24 (Test data)
0.69 (BCF) 3 (BCF; BCFWIN) -0.24 (Test data)
3 (BCF; BCFWIN) -0.24 (Test data)
3 (BCF; BCFWIN) -0.24 (Test data)
-0.24 (Test data)
` '
Not bioaccumulative.
501.187 (BCF; Other; Pimephales promelas)
3.5 - 3.94 (Calculated)
Potential for bioaccumulation (500 ≤ BCF ≤ 5000).
0.037 N/m (24 °C)
Toxic to bees. Not toxic to plants.
0.039 N/m (25 °C)
Toxic to bees.
0.0237 N/m
0.018 N/m (25 °C; 1 g/l)
Koc,2187.76; QSAR; log Koc; 3.34; QSAR
ment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Additional information : Handle empty containers with care because residual vapors are flammable.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

UN-No. (ADR) : 1993 : 1993 UN-No.(IATA)

14.2. UN proper shipping name

Proper Shipping Name (ADR) : FLAMMABLE LIQUID, N.O.S. Proper Shipping Name (IATA) : FLAMMABLE LIQUID, N.O.S.

Transport document description (ADR) UN 1993 FLAMMABLE LIQUID, N.O.S. (FLAMMABLE LIQUID, N.O.S.), 3, II, (D/E),

ENVIRONMENTALLY HAZARDOUS

14.3. Packing group

Class (ADR) : 3 Classification code (ADR) : F1 Class (IATA) : 3

29/02/2016 EN (English US) 7/9

Safety Data Sheet

Hazard labels (ADR) : 3



Hazard labels (IATA) : 3



14.4. Packing group

Packing group (ADR) : II Packing group (IATA) : II

14.5. Environmental hazards

Dangerous for the environment



Other information : No supplementary information available.

14.6. Special precautions for user

14.6.1. Overland transport

Hazard identification number (Kemler No.) : 33 Classification code (ADR) : F1

Orange plates :

33 1993

Special provision (ADR) : 274, 601, 640D

Transport category (ADR) : 2
Tunnel restriction code (ADR) : D/E
Limited quantities (ADR) : 11
Excepted quantities (ADR) : E2

14.6.2. Transport by sea

No additional information available

14.6.3. Air transport

CAO packing instructions (IATA) : 364 CAO max net quantity (IATA) : 60L PCA packing instructions (IATA) : 353 PCA Limited quantities (IATA) : Y341 PCA limited quantity max net quantity (IATA) : 1L PCA max net quantity (IATA) : 5L PCA Excepted quantities (IATA) : E2 ERG code (IATA) : 3H

14.6.4. Inland waterway transport

Carriage prohibited (ADN) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

29/02/2016 EN (English US) 8/9

Safety Data Sheet

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending

Regulation (EC) No 1907/2006.

Other information : None.

PHV SDS EU

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.

29/02/2016 EN (English US) 9/9