

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Date of issue: 18/04/2018

Revision date:

Version: 1.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Product name : 1,3-Dinitrobenzene SS
Product code : AL0-130286
Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Laboratory Use
Industrial/Professional use spec : For professional use only
Industrial
Use of the substance/mixture : Certified reference material for laboratory use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Phenova
6390 Joyce Dr. Suite 100
80403 Golden, CO - United States
T 1-866-942-2978 - F 1-866-283-0269
info@phenova.com - www.phenova.com

1.4. Emergency telephone number

Emergency number : ChemTel Assistance (US/Canada) 1-800-255-3924
ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute Tox. 4 (Oral) H302
Acute Tox. 4 (Dermal) H312
Acute Tox. 4 (Inhalation) H332
Carc. 2 H351

Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]

Carc.Cat.3; R40
Xn; R20/21/22
Full text of R-phrases: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS07

GHS08

Signal word (CLP) : Warning
Hazardous ingredients : 1,3-dinitrobenzene; Methylene Chloride
Hazard statements (CLP) : H302+H312+H332 - Harmful if swallowed, in contact with skin or if inhaled
H351 - Suspected of causing cancer

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Precautionary statements (CLP) : P261 - Avoid breathing dust/fume/gas/mist/vapors/spray
P301+P330+P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting
P270 - Do not eat, drink or smoke when using this product
P271 - Use only outdoors or in a well-ventilated area
P280 - Wear protective gloves, protective clothing, eye protection, face protection
P303+P361+P353+P310 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor
P304+P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing
P308+P313 - IF exposed or concerned: Get medical advice/attention
P362+P364 - Take off contaminated clothing and wash it before reuse
P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

No labeling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/Information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

| Name | Product identifier | % | Classification according to Regulation (EC) No. 1272/2008 [CLP] |
|--------------------|--|------|--|
| Methylene Chloride | (CAS No) 75-09-2 (EC-No.) 200-838-9 (EC index no) 602-004-00-3 | 99.8 | Carc. 2, H351 |
| 1,3-dinitrobenzene | (CAS No) 99-65-0 (EC-No.) 202-776-8 (EC index no) 609-004-00-2 | 0.2 | Acute Tox. 2 (Oral), H300 Acute Tox. 1 (Dermal), H310 Acute Tox. 2 (Inhalation), H330 STOT RE 2, H373 Aquatic Acute 1, H400 Aquatic Chronic 1, H410 |

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. IF exposed or concerned: Get medical advice/attention.

First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.

First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse.

First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persists.

First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

No additional information available

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.

Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect spillage. Store away from other materials.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood.

Hygiene measures : Gently wash with plenty of soap and water. Remove/Take off immediately all contaminated clothing. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container closed when not in use. Keep container tightly closed and in a well-ventilated place. Keep away from any flames or sparking source.

Incompatible products : Strong bases. Strong acids.

Incompatible materials : Sources of ignition. Direct sunlight.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

| Methylene Chloride (75-09-2) | | |
|------------------------------|----------------------------------|---|
| Belgium | Limit value (mg/m ³) | 177 mg/m ³ (Chlorure de méthylène; Belgium; Time-weighted average exposure limit 8 h) |
| Belgium | Limit value (ppm) | 50 ppm (Chlorure de méthylène; Belgium; Time-weighted average exposure limit 8 h) |
| France | VLE (mg/m ³) | 356 mg/m ³ (Dichlorométhane; France; Short time value; VRC: Valeur réglementaire contraignante) |
| France | VLE (ppm) | 100 ppm (Dichlorométhane; France; Short time value; VRC: Valeur réglementaire contraignante) |
| France | VME (mg/m ³) | 178 mg/m ³ (Dichlorométhane; France; Time-weighted average exposure limit 8 h; VRC: Valeur réglementaire contraignante) |
| France | VME (ppm) | 50 ppm (Dichlorométhane; France; Time-weighted average exposure limit 8 h; VRC: Valeur réglementaire contraignante) |
| Italy - Portugal - USA ACGIH | ACGIH TWA (ppm) | 50 ppm (Dichloromethane (Methylene chloride); USA; Time-weighted average exposure limit 8 h; TLV - Adopted Value) |
| United Kingdom | WEL TWA (mg/m ³) | 350 mg/m ³ Dichloromethane; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL TWA (ppm) | 100 ppm Dichloromethane; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL STEL (mg/m ³) | 1060 mg/m ³ Dichloromethane; United Kingdom; Short time value; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL STEL (ppm) | 300 ppm Dichloromethane; United Kingdom; Short time value; Workplace exposure limit (EH40/2005) |

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

| 1,3-dinitrobenzene (99-65-0) | | |
|------------------------------|----------------------------------|---|
| Belgium | Limit value (mg/m ³) | 1 mg/m ³ (Dinitrobenzène (tous isomères); Belgium; Time-weighted average exposure limit 8 h) |
| Belgium | Limit value (ppm) | 0.15 ppm (Dinitrobenzène (tous isomères); Belgium; Time-weighted average exposure limit 8 h) |
| France | VME (mg/m ³) | 1 mg/m ³ (Dinitrobenzène (tous isomères); France; Time-weighted average exposure limit 8 h; VL: Valeur non réglementaire indicative) |
| France | VME (ppm) | 0.15 ppm (Dinitrobenzène (tous isomères); France; Time-weighted average exposure limit 8 h; VL: Valeur non réglementaire indicative) |
| Italy - Portugal - USA ACGIH | ACGIH TWA (ppm) | 0.15 ppm (Dinitrobenzene, all isomers; USA; Time-weighted average exposure limit 8 h; TLV - Adopted Value) |
| United Kingdom | WEL TWA (mg/m ³) | 1 mg/m ³ Dinitrobenzene, all isomers; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL TWA (ppm) | 0.15 ppm Dinitrobenzene, all isomers; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL STEL (mg/m ³) | 3.5 mg/m ³ Dinitrobenzene, all isomers; United Kingdom; Short time value; Workplace exposure limit (EH40/2005) |
| United Kingdom | WEL STEL (ppm) | 0.5 ppm Dinitrobenzene, all isomers; United Kingdom; Short time value; Workplace exposure limit (EH40/2005) |


8.2. Exposure controls

Appropriate engineering controls

: Either local exhaust or general room ventilation is usually required.

Personal protective equipment

: Avoid all unnecessary exposure. Gloves. Protective clothing. Protective goggles. Safety glasses.



Hand protection

: Wear chemically resistant protective gloves. Wear suitable gloves resistant to chemical penetration.

Eye protection

: Chemical goggles or safety glasses. Safety glasses.

Skin and body protection

: Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin contact.

Respiratory protection

: Wear appropriate mask.

Other information

: Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| | |
|---------------------------|---------------------|
| Physical state | : Liquid |
| Color | : Colorless. |
| Odor | : characteristic. |
| pH | : No data available |
| Melting point | : No data available |
| Freezing point | : No data available |
| Boiling point | : No data available |
| Flash point | : No data available |
| Auto-ignition temperature | : No data available |
| Decomposition temperature | : No data available |
| Flammability (solid, gas) | : Non flammable. |
| Relative density | : No data available |
| Solubility | : No data available |
| Explosive properties | : No data available |
| Oxidizing properties | : No data available |
| Explosion limits | : No data available |

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Not established.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

fume. Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Oral: Harmful if swallowed. Dermal: Harmful in contact with skin. Inhalation: Harmful if inhaled.

| 1,3-Dinitrobenzene SS | |
|------------------------------|---|
| ATE CLP (oral) | 500 mg/kg body weight |
| ATE CLP (dermal) | 1100 mg/kg body weight |
| ATE CLP (gases) | 4500 ppmV/4h |
| ATE CLP (vapors) | 11 mg/l/4h |
| ATE CLP (dust, mist) | 1.5 mg/l/4h |
| Methylene Chloride (75-09-2) | |
| LD50 oral rat | > 2000 mg/kg (Rat; Literature study) |
| LD50 dermal rabbit | > 2000 mg/kg (Rabbit; Literature study) |
| 1,3-dinitrobenzene (99-65-0) | |
| LD50 oral rat | 60 mg/kg (Rat) |
| LD50 dermal rat | 1200 mg/kg (Rat) |
| ATE CLP (oral) | 5 mg/kg body weight |
| ATE CLP (dermal) | 5 mg/kg body weight |
| ATE CLP (gases) | 100 ppmV/4h |
| ATE CLP (vapors) | 0.5 mg/l/4h |
| ATE CLP (dust, mist) | 0.05 mg/l/4h |

| | |
|--|--|
| Skin corrosion/irritation | : Not classified Based on available data, the classification criteria are not met |
| Serious eye damage/irritation | : Not classified Based on available data, the classification criteria are not met |
| Respiratory or skin sensitization | : Not classified Based on available data, the classification criteria are not met |
| Germ cell mutagenicity | : Not classified Based on available data, the classification criteria are not met |
| Carcinogenicity | : Suspected of causing cancer. May cause cancer |
| Reproductive toxicity | : Not classified Based on available data, the classification criteria are not met |
| Specific target organ toxicity – single exposure | : Not classified Based on available data, the classification criteria are not met |
| Specific target organ toxicity – repeated exposure | : Not classified Based on available data, the classification criteria are not met |
| Aspiration hazard | : Not classified Based on available data, the classification criteria are not met |

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Potential Adverse human health effects and symptoms : Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1. Toxicity

| Methylene Chloride (75-09-2) | |
|------------------------------|--|
| LC50 fish 1 | 193 mg/l (LC50; 96 h; Pimephales promelas) |
| EC50 Daphnia 1 | 168.2 mg/l (EC50; 48 h) |
| 1,3-dinitrobenzene (99-65-0) | |
| LC50 fish 1 | 1.7 mg/l (LC50; 96 h) |
| EC50 Daphnia 1 | 27.4 mg/l (EC50; 48 h) |

12.2. Persistence and degradability

| 1,3-Dinitrobenzene SS | |
|-------------------------------|---|
| Persistence and degradability | Not established. |
| Methylene Chloride (75-09-2) | |
| Persistence and degradability | Not readily biodegradable in water. Biodegradable in the soil. |
| 1,3-dinitrobenzene (99-65-0) | |
| Persistence and degradability | Not readily biodegradable in water. Non degradable in the soil. |

12.3. Bioaccumulative potential

| 1,3-Dinitrobenzene SS | |
|------------------------------|--|
| Bioaccumulative potential | Not established. |
| Methylene Chloride (75-09-2) | |
| BCF fish 1 | 2 - 40 (BCF) |
| Log Pow | 1.25 (Experimental value) |
| Bioaccumulative potential | Low potential for bioaccumulation (BCF < 500). |
| 1,3-dinitrobenzene (99-65-0) | |
| BCF fish 1 | 4.5 - 7.5 (BCF; 72 h) |
| BCF fish 2 | 74.13 (BCF) |
| Log Pow | 1.49 - 1.6 |
| Bioaccumulative potential | Low potential for bioaccumulation (BCF < 500). |

12.4. Mobility in soil

| Methylene Chloride (75-09-2) | |
|------------------------------|---|
| Surface tension | 0.028 N/m (20 °C) |
| Ecology - soil | May be harmful to plant growth, blooming and fruit formation. |

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

UN-No. (ADR) : 2810
UN-No. (IATA) : 2810
UN-No. (IMDG) : 2810
UN-No. (ADN) : 2810

14.2. UN proper shipping name

Proper Shipping Name (ADR) : TOXIC LIQUID, ORGANIC, N.O.S.
Proper Shipping Name (IATA) : Toxic liquid, organic, n.o.s.
Proper Shipping Name (IMDG) : TOXIC LIQUID, ORGANIC, N.O.S.
Proper Shipping Name (ADN) : TOXIC LIQUID, ORGANIC, N.O.S.

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Transport document description (ADR) : UN 2810 TOXIC LIQUID, ORGANIC, N.O.S., 6.1, III, (E)

14.3. Packing group

Class (ADR) : 6.1
Classification code (ADR) : T1
Class (IATA) : 6.1
Class (IMDG) : 6.1
Class (ADN) : 6.1
Classification code (ADN) : T1
Hazard labels (ADR) : 6.1



Division (IATA) : 6.1
Hazard labels (IATA) : 6.1



Hazard labels (IMDG) : 6.1



Hazard labels (ADN) : 6.1



14.4. Packing group

Packing group (ADR) : III
Packing group (IATA) : III
Packing group (IMDG) : III
Packing group (ADN) : III


14.5. Environmental hazards

Other information : No supplementary information available.

14.6. Special precautions for user

14.6.1. Overland transport

Hazard identification number (Kemler No.) : 60
Classification code (ADR) : T1
Orange plates :



Special provision (ADR) : 274, 614
Transport category (ADR) : 2
Tunnel restriction code (ADR) : E
Limited quantities (ADR) : 5I
Excepted quantities (ADR) : E1
EAC : 2X
APP : B

14.6.2. Transport by sea

Special provision (IMDG) : 223, 274
Limited quantities (IMDG) : 5 L
Excepted quantities (IMDG) : E1

1,3-Dinitrobenzene SS

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

| | |
|------------------------------------|---|
| Packing instructions (IMDG) | : P001, LP01 |
| IBC packing instructions (IMDG) | : IBC03 |
| Tank instructions (IMDG) | : T7 |
| Tank special provisions (IMDG) | : TP1, TP28 |
| EmS-No. (Fire) | : F-A |
| EmS-No. (Spillage) | : S-A |
| Stowage category (IMDG) | : A |
| Properties and observations (IMDG) | : Toxic if swallowed, by skin contact or by inhalation. |

14.6.3. Air transport

| | |
|--|----------------|
| CAO packing instructions (IATA) | : 663 |
| CAO max net quantity (IATA) | : 220L |
| PCA packing instructions (IATA) | : 655 |
| PCA Limited quantities (IATA) | : Y642 |
| PCA limited quantity max net quantity (IATA) | : 2L |
| PCA max net quantity (IATA) | : 60L |
| PCA Excepted quantities (IATA) | : E1 |
| Special provision (IATA) | : A3, A4, A137 |
| ERG code (IATA) | : 6L |

14.6.4. Inland waterway transport

| | |
|-----------------------------------|------------------|
| Special provision (ADN) | : 274, 614, 802 |
| Limited quantities (ADN) | : 5 L |
| Excepted quantities (ADN) | : E1 |
| Carriage permitted (ADN) | : T |
| Equipment required (ADN) | : PP, EP, TOX, A |
| Ventilation (ADN) | : VE02 |
| Number of blue cones/lights (ADN) | : 0 |
| Carriage prohibited (ADN) | : No |

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions
Contains no REACH candidate substance
Contains no REACH Annex XIV substances.

15.1.2. National regulations

Germany

Water hazard class (WGK) : 2 - hazardous to water

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

| | |
|-------------------|--|
| Data sources | : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. |
| Other information | : None. |

PHV SDS EU

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.