

Safety Data Sheet

according to Regulation (EC) No. 453/2010

Date of issue: 08/04/2014 Revision date: 14/04/2015 : Version: 1.1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : N-Nitrosodiphenylamine Standard

Product code : AL0-101245
Product group : Trade product

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Laboratory Use Industrial/Professional use spec : Industrial

For professional use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Phenova

6390 Joyce Dr. Suite 100 80403 Golden, CO - United States T 1-866-942-2978 - F 1-866-283-0269 info@phenova.com - www.phenova.com

1.4. Emergency telephone number

Emergency number : ChemTel Assistance (US/Canada) 1-800-255-3924 ChemTel Assistance (International) +1 813-248-0585

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Carc. 2 H351

Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]

Carc.Cat.3; R40

Full text of R-phrases: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labeling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :

GHS08

Signal word (CLP) : Warning

Hazardous ingredients : Methylene Chloride

Hazard statements (CLP) : H351 - Suspected of causing cancer

Precautionary statements (CLP) : P281 - Use personal protective equipment as required

P308+P313 - IF exposed or concerned: Get medical advice/attention

P405 - Store locked up

No labeling applicable

2.3. Other hazards

No additional information available

16/04/2015 EN (English US) 1/1

Safety Data Sheet

according to Regulation (EC) No. 453/2010

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Methylene Chloride (Component)	(CAS No) 75-09-2 (EC no) 200-838-9 (EC index no) 602-004-00-3	99.9	Carc. 2, H351
N-nitrosodiphenylamine (Component)	(CAS No) 86-30-6 (EC no) 201-663-0	0.1	Acute Tox. 4 (Oral), H302 Carc. 2, H351 Aquatic Chronic 2, H411

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. IF exposed or concerned: Get

medical advice/attention.

First-aid measures after inhalation : Allow victim to breathe fresh air. Allow the victim to rest.

First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed

by warm water rinse.

First-aid measures after eye contact : Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness

persist.

First-aid measures after ingestion : Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

No additional information available

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.

Emergency procedures : Ventilate area.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect

spillage. Store away from other materials.

6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

16/04/2015 EN (English US) 2/1

Safety Data Sheet

according to Regulation (EC) No. 453/2010

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or

smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. Obtain special instructions before use. Do not handle until all safety precautions have

been read and understood.

Hygiene measures : Gently wash with plenty of soap and water. Remove/Take off immediately all contaminated

clothing. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container closed when not in use. Keep container tightly closed and in a well-ventilated

place. Keep away from any flames or sparking source.

Incompatible products : Strong bases. Strong acids.
Incompatible materials : Sources of ignition. Direct sunlight.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Appropriate engineering controls : Either local exhaust or general room ventilation is usually required.

Personal protective equipment : Avoid all unnecessary exposure. Gloves. Protective clothing. Protective goggles. Safety

glasses.







Hand protection : Wear chemically resistant protective gloves. Wear suitable gloves resistant to chemical

penetration.

Eye protection : Chemical goggles or safety glasses. Safety glasses.

Skin and body protection : Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin

contact.

Respiratory protection : Wear appropriate mask.

Other information : Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid Color Colorless Odor characteristic. рΗ No data available Melting point : No data available Freezing point No data available **Boiling point** No data available Flash point No data available Auto-ignition temperature No data available No data available Decomposition temperature Flammability (solid, gas) Non flammable Relative density No data available Solubility No data available Explosive properties No data available Oxidizing properties No data available : No data available **Explosion limits**

9.2. Other information

No additional information available

16/04/2015 EN (English US) 3/1

Safety Data Sheet

according to Regulation (EC) No. 453/2010

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Not established.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

fume. Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

N-nitrosodiphenylamine (86-30-6)		
LD50 oral rat	1650 mg/kg (Rat)	
LD50 dermal rabbit	> 7940 mg/kg (Rabbit)	
ATE CLP (oral)	1650.000 mg/kg body weight	

Methylene Chloride (75-09-2)	
LD50 oral rat	> 2000 mg/kg (Rat; Literature study)
LD50 dermal rabbit	> 2000 mg/kg (Rabbit; Literature study)

Skin corrosion/irritation : Not classified

Based on available data, the classification criteria are not met

Serious eye damage/irritation : Not classified

Based on available data, the classification criteria are not met

Respiratory or skin sensitization : Not classified

Based on available data, the classification criteria are not met

Germ cell mutagenicity : Not classified

Based on available data, the classification criteria are not met

Carcinogenicity : Suspected of causing cancer.

May cause cancer

Reproductive toxicity : Not classified

Based on available data, the classification criteria are not met

Specific target organ toxicity (single exposure) : Not classified

Based on available data, the classification criteria are not met

Specific target organ toxicity (repeated

exposure)

: Not classified

Based on available data, the classification criteria are not met

Aspiration hazard : Not classified

Based on available data, the classification criteria are not met

Potential Adverse human health effects and

symptoms

Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1. Toxicity

N-nitrosodiphenylamine (86-30-6)		
LC50 fish 1	6.4 mg/l (48 h; Oryzias latipes; Static system)	
EC50 Daphnia 1	7.8 mg/l (48 h; Daphnia magna)	
LC50 fish 2	5.8 mg/l (96 h; Lepomis macrochirus)	
Methylene Chloride (75-09-2)		
LC50 fish 1	193 mg/l (96 h; Pimephales promelas; Flow-through system)	
EC50 Daphnia 1	168.2 mg/l (48 h; Daphnia magna)	
LC50 fish 2	220 mg/l (96 h; Lepomis macrochirus; Flow-through system)	
Threshold limit algae 1	1450 mg/l (192 h: Scenedesmus quadricauda: Cell numbers)	

16/04/2015 EN (English US) 4/1

N-Nitrosodiphenylamine Standard Safety Data Sheet according to Regulation (EC) No. 453/2010

Class (IMDG)

Class (ADN)

W (I I OII II (77 00 0)	
Methylene Chloride (75-09-2)	770 # (400 L M)
Threshold limit algae 2	550 mg/l (192 h; Microcystis aeruginosa)
12.2. Persistence and degradability	
N-Nitrosodiphenylamine Standard	
Persistence and degradability	Not established.
·	1100 COLUMNOTICU.
N-nitrosodiphenylamine (86-30-6)	Notice district and delicity with Disdensed billion in 20 and delicity
Persistence and degradability	Not readily biodegradable in water. Biodegradability in soil: no data available.
Methylene Chloride (75-09-2)	
Persistence and degradability	Not readily biodegradable in water. Biodegradable in the soil.
12.3. Bioaccumulative potential	
N-Nitrosodiphenylamine Standard	
Bioaccumulative potential	Not established.
N-nitrosodiphenylamine (86-30-6)	
BCF fish 1	217 (336 h; Lepomis macrochirus; Flow-through system)
BCF fish 2	4.6 - 38 (Cyprinus carpio; Test duration: 6 weeks)
Log Pow	3.13 - 3.96
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).
Methylene Chloride (75-09-2)	
BCF fish 1	2 - 40 (Cyprinus carpio; Test duration: 6 weeks)
Log Pow	1.25 (Experimental value)
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).
12.4. Mobility in soil	
Methylene Chloride (75-09-2)	0.000 N/m (00 °C)
Surface tension	0.028 N/m (20 °C)
Ecology - soil	May be harmful to plant growth, blooming and fruit formation.
12.5. Results of PBT and vPvB assessme	ent en
No additional information available	
12.6. Other adverse effects	
Additional information	: Avoid release to the environment
SECTION 42: Discussed associations tile	**
SECTION 13: Disposal consideratio	ns
13.1. Waste treatment methods	
Waste disposal recommendations	: Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials	: Avoid release to the environment.
SECTION 14: Transport information	
In accordance with ADR / RID / IMDG / IATA / A	
14.1. UN number	
UN-No. (ADR)	: 2810
, ,	: 2810 : 2810
UN-No.(IATA)	. 2010
14.2. UN proper shipping name	
Proper Shipping Name (ADR)	: TOXIC LIQUID, ORGANIC, N.O.S.
Proper Shipping Name (IATA)	: TOXIC LIQUID, ORGANIC, N.O.S.
Proper Shipping Name (IMDG)	: TOXIC LIQUID, ORGANIC, N.O.S.
Proper Shipping Name (ADN)	: TOXIC LIQUID, ORGANIC, N.O.S.
Transport document description (ADR)	: UN 2810 TOXIC LIQUID, ORGANIC, N.O.S. (dichloromethane(75-09-2)), 6.1, III, (E)
14.3. Packing group	
Class (ADR)	: 6.1
Classification code (ADR)	: T1
, ,	: 6.1
Class (IATA)	. 0.1

16/04/2015 EN (English US) 5/1

: 6.1

: 6.1

Safety Data Sheet

according to Regulation (EC) No. 453/2010

Hazard labels (ADR) : 6.1

6

Hazard labels (IATA) : 6.1



14.4. Packing group

Packing group (ADR) : III
Packing group (IATA) : III

14.5. Environmental hazards

Other information : No supplementary information available.

14.6. Special precautions for user

14.6.1. Overland transport

Hazard identification number (Kemler No.) : 60
Classification code (ADR) : T1

Orange plates : 6

60 2810

Special provision (ADR) : 274, 614

Transport category (ADR) : 2
Tunnel restriction code (ADR) : E
Limited quantities (ADR) : 51
Excepted quantities (ADR) : E1

14.6.2. Transport by sea

No additional information available

14.6.3. Air transport

CAO packing instructions (IATA) : 663 CAO max net quantity (IATA) : 220L : 655 PCA packing instructions (IATA) PCA Limited quantities (IATA) : Y642 PCA limited quantity max net quantity (IATA) : 2L PCA max net quantity (IATA) : 60L PCA Excepted quantities (IATA) : E1 Special provision (IATA) : A137 ERG code (IATA) : 6L

14.6.4. Inland waterway transport

Carriage prohibited (ADN) : No

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions

Contains no REACH candidate substance

Contains no REACH Annex XIV substances.

15.1.2. National regulations

No additional information available

16/04/2015 EN (English US) 6/1

Safety Data Sheet

according to Regulation (EC) No. 453/2010

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending

Regulation (EC) No 1907/2006.

Other information : None.

PHV SDS EU

Copyright 2015 Phenova, Inc. License granted to make paper copies for internal use. The information contained in this Safety Data Sheet is based on our current knowledge. The information contained in this document should be used only as a guide for appropriate safety precautions and should not be considered to be all inclusive. Users should make their own investigation to determine the suitability of the information for their particular purposes. The document does not represent any guarantee of the properties of the product. Phenova, Inc. shall not be held liable for any damage resulting from the handling or use of this product. Visit the Terms and Conditions of Sale link at www.phenova.com for additional terms and conditions of sale.

16/04/2015 EN (English US) 7/1